

## Product datasheet for **MC205753**

### Dalrd3 (NM\_026378) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Dalrd3 (NM\_026378) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Dalrd3  
**Synonyms:** 6330580J24Rik; C77829; mir-425  
**Mammalian Cell Selection:** Neomycin  
**Vector:** PCMV6-Kan/Neo (PCMV6KN)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Fully Sequenced ORF:** >BC013483

```

ATGGCGACCGGGCGCCTTGGGGTTGGCGAGACGCTGGAGGCCCTTAATGCAGCTGTAGGGCCTGGTAGCC
CAGTGTGGTTTAAAGGAAACCCATGCCCGGCACCTTCGTGTTTCGAGACTTCTTGGCACCGGAAGCGCGCT
ACAGGCGCGCTTACAGGGACGGCCAGGTGCCAGAGTGCCTGTTCCGTGCGGTCTCCTGCCTGCAGGGTCCC
GGAGTGGCCCTGTACTGCGCTGCGCGCCGACACCCGCGGGTTTATCGCTCCAACCTACAGCGGCCCGCTG
TCTTCAACATGTACTCGGTGCCCTGGCCTCCTATGCTACTCCCGCAAACCTGCCTCACCAGGCCCGCG
AGTTGTCTACACTGCCAGCACTGCGCTGTAATTCGGGGACGCTCCGATTGAGCCAGCTGCGTGCCGTG
CTCGTAGCCGATCATCTAGTACGGTACTGCGTGTCTATGGGGTACGTGTGTGCTCAGTGCCCCCGGTGC
GGGACCCTCACATGTCAACTTTCCTCCAGAACTTCGGGTGGAATGGCCACTGCCTCGAAGAGCACCTC
GACCGAAACACTGAGAACCTGTGTTCTTGCAAACCTTAATGGTCTAAGGAAGCGACGCTGCCCCCTGGC
GTCCTAGGCCGACTGTGCTGAAGGAACGGTAGAACAGAGAGGCACTGCTGGCTATGACCCAGCATAG
ACCACTGTTTGGTGAAGATGTGCTCTCTGTGCTCTCGGAGCTGCAGGAGGCTGTGCCCACTGGCC
AGAGGGCGGCCACCCAGGCCCGGTGGGCGTCCAGATGCTGGTGTAGATGACTGCGTGGTTATACACGTG
GTTAGCTGTGAAGAGGCATTCCAGCAACAAAAGCTGGATTTGCTCTGGCAGAAATTGGATGACCGGGCTC
CTCACAAGCAGAAACACCTGGTGTGTGGTCCCCTGAAAATGGCTGGTGTACCCGGCACCCAGCTGACTGC
CCCTCAGTACTACAGGCTCAGACATGCTCAGGTGTGCGAGGCCTCAGCTCTGAAGCACGGAGGAGACCTG
GCACAAGACCCCGCTGGACTGAGACCTTCGACATCCTCTCCGTGGTACCATCAAATTTGAGATGCTAA
GCACAGCCCCACAGAGCCAGCTCCTCCTGGCCACAGCACCATCTCAACAAGGGCACAAGAGCGGGAC
TTTTGTATGTATAACTGTGCCCGCCTTGCCACACTCTTCGAGGGCTACAACATGGAACGGAGCAAGGT
CTGTATCCCACTTTCCGCTTGTGAGCAGTCTGGACTTCTCACTACTCCATGATGAGGGTGAAGTGGCTGC
TGCTTTTCAACAGTGTCTTCCCTTTCTGGACCTGCTGAGCCAGACTGTGAGTCTGGCCGACCCCAAGG
CCTCCACATCCCTGTCCGCACGGAGATGGTATGCAAGTTCCTGGTGCAGCTCAGCATGGATTTACGCTCC
TACTATAACCGAGTACATATTCTAGGGAACTCGGCCACACCTCTTTGGTCAAGTGTGCTCGCTAC
AGCTTTTGAAGGCGAGTGCCTGAGGTGTTCCACACAGGCTTGCTATGCTGGGCTCCCTCCGCTGAGCCA
CATCTAAGGGCCAATAACTGGGAATGTTATAAAAAGCATCAACTGGGAAAGAAACAAGACCTCACAGAT
GTTGGCACTTAAAGTCAGAATAAATACAGTTTCATCCTATAAAAAAAAAAAAAAAAAA

```



[View online >](#)

<b>Restriction Sites:</b>	RsrII-NotI
<b>ACCN:</b>	NM_026378
<b>Insert Size:</b>	1617 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">BC013483</a> , <a href="#">AAH13483</a>
<b>RefSeq Size:</b>	1736 bp
<b>RefSeq ORF:</b>	1617 bp
<b>Locus ID:</b>	67789
<b>UniProt ID:</b>	<a href="#">Q6PJN8</a>
<b>Cytogenetics:</b>	9 F2
<b>Gene Summary:</b>	Involved in tRNA methylation. Facilitates the recognition and targeting of tRNA(Arg)(CCU) and tRNA(Arg)(UCU) substrates for N(3)-methylcytidine modification by METTL2.[UniProtKB/Swiss-Prot Function]