

Product datasheet for **MC205754**

Fah (NM_010176) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Fah (NM_010176) Mouse Untagged Clone
Tag: Tag Free
Symbol: Fah
Synonyms: swst
Mammalian Cell Selection: Neomycin
Vector: PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >BC010767
 CAGCATGTCCTTTATTCCAGTGGCCGAGGACTCCGACTTTCCCATCCAAAACCTGCCTATGGTGTITTC
 TCCACTCAAAGCAACCCAAAGCCACGGATTGGTGTAGCCATCGGTGACCAGATCTTGACCTGAGTGTCA
 TTAACACCTCTTTACCGGACCTGCCCTTTCCAAACATCAACATGTCTTCGATGAGACAACCTCTCAATAA
 CTTTCATGGGTCTGGGTCAAGCTGCATGGAAGGAGGCAAGAGCATCCTTACAGAACTTACTGTCTGCCAGC
 CAAGCCCGGCTCAGAGATGACAAGGAGCTTCGGCAGCGTGCATTACCTCCCAGGCTTCTGCGACAATGC
 ACCTTCTGCTACCATAGGAGACTACACGGACTTCTACTCTTCTCGGCAGCATGCCACCAATGTTGGCAT
 TATGTTTCAGAGGCAAGGAGAATGCGCTGTTGCCAAATTGGCTCCACTTACCTGTGGGATACCATGGCCGA
 GCTTCTCCATTGGTATCTGGAACCCCGATTGGAAGACCCATGGGGCAGATGAGACCTGATAACTCAA
 AGCCTCTGTGTATGGTGCCTGCAGACTCTTAGACATGGAGTTGGAATGGCTTTCTTCGTAGGCCCTGG
 GAACAGATTCGGAGAGCCAATCCCATTTCCAAAGCCCATGAACACATTTTCGGGATGGTCTCATGAAC
 GACTGGAGCGCACGAGACATCCAGCAATGGGAGTACGTCCCACTTGGGCCATTCTGGGGAAAAGCTTTG
 GAACCACAATCTCCCCGTGGGTGGTGCCTATGGATGCCCTCATGCCCTTTGTGGTGCCAAAACCCAAAGCA
 GGACCCCAAGCCCTTGCCATATCTTGCCACAGCCAGCCCTACACATTTGATATCAACCTGTCTGTCTCT
 TTGAAAGGAGAAGGAATGAGCCAGGCGGCTACCATCTGCAGGTCTAACTTTAAGCACATGTACTGGACCA
 TGCTGCAGCAACTCACACACCACTCTGTTAATGGATGCAACCTGAGACCTGGGGACCTCTTGCTTCTGG
 AACCATCAGTGGATCAGACCCTGAAAGCTTTGGCTCCATGCTGGAAGTGTCTGGAAGGGGAACAAAGGCC
 ATCGATGTGGGGCAGGGCAGACCAGGACCTTCTGCTGGACGGCGATGAAGTCATATAACAGGTCACCT
 GCCAGGGGACGGCTACCGTGTGGCTTTGGCCAGTGTGCTGGGAAAGTGTGCTGCCCTTTCCACCAGC
 CTGAAGCTCCGGAAGTCAAGACACACCTTGCCTTATGAGGATCATGCTACCCTGCATCAGTCAGGA
 ATGAATAAAGCTACTTTGATTGTGGGAAAAAAAAAAAAAAAAAAAA

Restriction Sites: RsrII-NotI
ACCN: NM_010176
Insert Size: 1260 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>BC010767</u> , <u>AAH10767</u>
RefSeq Size:	1373 bp
RefSeq ORF:	1260 bp
Locus ID:	14085
UniProt ID:	<u>P35505</u>
Cytogenetics:	7 48.36 cM