

Product datasheet for MC205767

Rbm4b (NM_025717) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Rbm4b (NM_025717) Mouse Untagged Clone
Tag: Tag Free
Symbol: Rbm4b
Synonyms: 4921506I22Rik; AI504630; AI506404; Lark2; Neosin
Mammalian Cell Selection: Neomycin
Vector: PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >BC019488

```

CTCTGACGCGGGACCCGCCGCCACCCGCCGCCACCCGGAGGCTCTTGTGAGGATGGTGAAGCTGTTCA
TTGGAAATCTGCCCGGGAGGCCACAGAGCAGGAGATCCGCTCACTCTTCGAGCAGTACGGGAAGGTGCT
GGAATGTGACATCATTAAGAATTATGGCTTTGTGCACATAGAGGACAAGACGGCCGCTGAGGATGCCATA
CGCAACCTGCACCACTACAAGCTGCACGGAGTGAACATCAATGTGGAAGCCAGCAAGAATAAGAGCAAAG
CTTCAACCAAGTTACACGTGGGCAACATCAGCCCCACTTGTACCAACCAAGAGCTTCGGGCCAAGTTTGA
GGAGTACGGCCAGTCATCGAATGTGACATCGTAAAAGATTATGCCTTTGTACACATGGAGCGGGCAGAG
GATGCGGTGGAGGCCATCAGGGGCTCGACAACACAGAGTTTCAAGGCAAACGAATGCATGTACAGTTGT
CAACTAGCCGGCTTCGGACTGCCCTGGGATGGGAGACCAGAGTGGTTGCTATCGGTGTGGTAAAGAAGG
ACATTGGTCCAAAGAGTCCCGTAGACCGTACAGGGCGTGTGGCAGACTTTACTGAGCAGTACAATGAA
CAATATGGAGCAGTGCACACACCTTATACTATGGGTTATGGGGAGTCCATGTATTACAACGATGCCTATG
GAGCACTTGACTACTATAAGCGCTACCGGTCCGATCTTACGAAGCTGTGGCGGGCGGCTGCAGCTTC
CGCATACAATTACGCAGAGCAGACCATGTCTCATCTTCTCAAGTTCAGAGCTCAGCTGTACCCAGTCAC
CTCAACTCCACTTCTGTTGATCCTTATGACAGACACCTATTGCAGAACTCTGGCTCTGCTGCCACCTCAG
CTGCAATGGCTGCTGCCGCTTCTCTTCTTATTATGGAAGGGACAGGAGCCCACTACGTCGGAATGCGGC
TGTGCTTCCCGCAGTTGGAGAGGGCTATGGTTATGGGCCAGAGAGCGAGATGTCTCAGGCTTCAGCAGCA
ACTCGGAATCTCTGTATGACATGGCCCGGTATGAGCGGGAGCAGTATGTGGACCGAACACGGTACTCAG
CCTTTTAAAACTGGAGGTAGGATAATTGCGGACTGAACCCTCGGGCTCGGTCATATATGAACTTGG
TCCTCGCGGTCCCCTTTGCCAGGATGTTCCATTGCTTCATGTTTCAGTAAACAAAGGAATTTGTACCA
ACTATGTTTTCTTCTTAATTTAATTCTTCTAAGTTGACTTTTCTTCTCGATGCTAGTTGTCTGTAGC
TTTTCACTGTTCTTATACCCTCAGCCTCTGAACAGCCCTAGGTAAGGGTTATGCTGACATCCCTTTTCC
TGTACAGTAGAAGCCCTCTTAATCTTGCTTTTCTTAGGAGTTGAGCCCTTCTCCCTGCCTTCTGCAGC
ATCTCCTTTCCCTTTAAATGACCATGTAGTGGCAAGCAACCTTAACTCTTCTGTTAGTGTGACTCT
TAGCATTGAAGCTAGTCTTCTGAAGTTGCTAGGACCATGGGTTTTGTTGTTGTTTTGTTTTGATTTTGT
TTTGTGTTTTCGGTTTTGCTGACCTGTGATCGTGGTACAGCATTTGCTGAAATTTAGCCTGTTTTATTCC
CACTCCTCCCAATTTTTTTTTGAAAAAAAAAATAAATGTTTCTAATACTTAAAAAAAAAAAAAAAAA

```



[View online »](#)

Restriction Sites:	RsrII-NotI
ACCN:	NM_025717
Insert Size:	1074 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC019488 , AAH19488
RefSeq Size:	1746 bp
RefSeq ORF:	1074 bp
Locus ID:	66704
UniProt ID:	Q8VE92
Cytogenetics:	19 A
Gene Summary:	Required for the translational activation of PER1 mRNA in response to circadian clock. Binds directly to the 3' UTR of the PER1 mRNA.[UniProtKB/Swiss-Prot Function]