

Product datasheet for MC205771

Syncn (NM_026716) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Syncn (NM_026716) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Syncn
Synonyms:	0910001K16Rik; 1810038B08Rik; SYL
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC019567 AGACCTCTCGTGGTTCATTGTCACCATGTCCCTGCTGTGCCCACTGCTGTGGCTTTAGCCCTGGTGGCT GTCCCTGGTGCCCATGGCAATTGCCCGGTGCCCGCAGACCTGAAGAAGTCAGATGGGACGCGTACGTGCG CCAGACTCTATGAGAAGAGTGATCCATACTATGACAACCTGCTGCCAGGGGCCGTTGTCTGTGGAGCC AGGTAAGTACCTGCCCTATCTGCCCTCGGGCTGGTCTAACTGTCATCGTCTCTGGTAGTGGGCCAGCGC TGTGAGATCACTGTGTGGTCTCTCCCGGTAAGCATGGCAAGACCGCAAGTTCACCGCGGGTACTTACC CTCGCCTGGAAGAGTACCGCAAAGGCATCTTTGGAGACTGGTCCGACTCCATCTCTGCCCTCTACTGCAA GTGCTATTGATGCCTTGAGGCTCTGTCTACCCAGCCTGGCCTTGGGAATTGCTGTAGCTCCAAGAGCCAG GAGGCAAGATGACCCACGACCTGCTCTCATAGCTTCCCTGTAATACAGCCCTTTCAAAGTTAAAAAAA AAAAAAAAAAAAAAAAAAAAAAAAAAAAAA
Restriction Sites:	RsrII-NotI
ACCN:	NM_026716
Insert Size:	405 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).


[View online »](#)

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>BC019567</u> , <u>AAH19567</u>
RefSeq Size:	589 bp
RefSeq ORF:	405 bp
Locus ID:	68416
UniProt ID:	<u>Q8VCK7</u>
Cytogenetics:	7 B1
Gene Summary:	Functions in exocytosis in pancreatic acinar cells regulating the fusion of zymogen granules with each other. May have a pore-forming activity on membranes and regulate exocytosis in other exocrine tissues.[UniProtKB/Swiss-Prot Function]