

Product datasheet for MC205823

Lzic (NM_026963) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Lzic (NM_026963) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Lzic
Synonyms:	1810030J04Rik; AW047580
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC023174

```

GGGGAGGGCGAAGGAGGTTGTCCTGTTTCCGATCAGGGAGGACGGAGAAGCGTTTCGCAGGCGTAGGACT
CCCCGCCTCCTCAGGTGACCACCGGGCCTGCTGGAGCCCGCCTTGGCGGACCTGCACACCTCGCCGGG
AGCTTGGAGTTCTGTCTGTCCGCGGGAGGAGCCGGCATCCGAGCGTGCGGTCCGGCTCTGCCTATTTT
TGCTCCGGTGGTCCCGGGCCCGGCTTTGCTCACCTCCGTCCTGAGAAGGTGTGCTTTCTGCCCTTGCC
TTCTCCTGCGTCTGAGTTGGACAGTTTTCTCAAGCCGAGTCCCGAGTGTGTGGCTCACCCAGAGGAG
ACTCGGCAAAAAGCTACGGGAGGCACAGCCAGGGATTCAAGTGGCCTCTCCCGCTTCTCAGGAATGGA
GAACGTATTGCAGTTACAGTCTGCTTAGCGGGAAGTGATCTGGACAGCTGTAGGCCAAGTTCAGACACAG
CAAATAAACTCTTAACAGTTTTACCTTGTCTAGCCAGTCCACCCAGGAGTTTGAAGTTTGAGCCT
CCGGTAACTTTATGACAACTTCGCTTTCTCAATTTGATGACACTCAGAGCTCCGACAGTTGGCGAGTT
TTGAGCACCGTGCACGGGACAGGACACTGGAAGCCGTGATTAATGGCTTCCAGAGGAAAGACAGAGA
CAAGCAAATAAAGCAGAATTTAGAAGAGCAGTTGGATAGACTGATGCAGCAGTTACAAGATCTGGAGGA
ATGCAGGGAGGAGCTCGATGCGGATGAATATGAAGAACTAAAAAGGAACTTTGGAGCAGCTGAGTGAG
TTCAACGATTTCGCTGAAGAAGATTATGTGAGGAAATATGACGTTGGTAGATGAACTCAGTGGGATGCAAC
TGGCCATCCAGGCAGCCATCAGCCAGGCCTTTAAACTCCTGAAGTATCAGATTGTTTGGCAAGAAACA
ACCAGGCAACTTCGGACAAGGCTAGCCGAGATGGACAGAGATCTGATGGTGGGGAAGCTGGAAAGAGAG
CTGTATACCCAGCGGAAGGTGGAGATACTGACAGCACTCAGGAACTCGGAGAGAAGCTGACAGAAGACG
ATGAGACCTTCTGTGTCAGCCAATGCAAGTGTGTGCTCAGCCAGTTTGAAGAACTCCACAGAGCTTGG
CTCTGGAGACAAGTCTTGTCTGCTGGCAGGGTTCGATGTTGAAAAAGCCAAGAAATGACGGCGTGGTGGG
TTTGGCTCTGCCGCCCCACCCCTCCTGCACACGGTGTCCATCCTGAGGACTCAGGCCACACAAGCT
GTAGAGACTCTCACAATGCAGCAACAGCCAGGAACAGGTGACTTGAACACAGTGTGGCTTCTGCCTAG
GGTCTCGTCCCTCTCTGTGGATCTCCCTTCTTACGTGACCCAGGCTTGTGGACAGAGGGACATTG
TTCACTATGATTCGTTTTGTTAAGATATCTTGGTTTTAGTGTGATAGGTGTGAGAGCTATGAGCTAAGA
CTGTCTGACCTCTGTCTGCATGCCTTGGGATGGACAGCAGTGTGTCCCTGGCCCTCCTGTGGCTGGAG
AATCTTTTCAGTTTTTATAACAAGGAAACCAAGACTTAGCACTTGTGATGTCTAAAGTCACATTAGAGGCA
GAGCTGAGGCTGGGGCTCAGGTGTCTGAATCCAAGTTTTACTTCTTTATCTGTACTGTTTTATATATA
TATAATAATTAAGTATTGTTGTGCAAAAAAAAAAAAAAAAAA

```



[View online »](#)

Restriction Sites:	RsrII-NotI
ACCN:	NM_026963
Insert Size:	573 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>BC023174</u> , <u>AAH23174</u>
RefSeq Size:	1791 bp
RefSeq ORF:	573 bp
Locus ID:	69151
UniProt ID:	<u>Q8K3C3</u>
Cytogenetics:	4 E2