

## Product datasheet for MC205839

### Hhatl (NM\_029095) Mouse Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Hhatl (NM_029095) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Hhatl
Synonyms:	1110011D13Rik; Gup1; Mg56
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC024464

```
AAGCCCCACAGCTCAGCCAGCCGGACAGACTCACGGTCGGACTCAAGACCCCGGAGCCCTGAGGTGGGCA
GGCGCCAGGGTTCCTCGCAGCCTCTTCAAGCCCTCCTCAGGCCATGGGCATCAAGACAGCACTGCCTGC
AGCTGAGCTCGGCCTCTATTCACCTGGTGTGAGCGGTGCCCTGGCTTATGCCGGCCGAGGTCTCCTGGAG
GCATCACAAAGATGGGGCCACAGGAAGGCCCTCCGGGAGTCTGTGAGACCTGGCTGGGAGTACCTGGGCC
GGAAGATGGATGTGGCTGACTTCGAGTGGGTGATGTGGTTTACCAACTCCGAAATGTCATTGTCTTTGC
CCTCTCTGGACACGTGCTGTTTGCCAACTCTGCACGATGGTGGCCCCCAGCTCCGCTCCTGGATGTAT
GCTGTATATGGAGTCTGGCTGTGGTGGGCACAATGGGCCATGGTACCTGCTGCTGCTGCTTGGCCACT
GTATGGTCTCTATGTGGCCTCACTCTTGGCCAGCGTTGGCTGTGCCTTGCCTTGGTCTGGCCAGCTT
GGCGTCTTCAAGTTGACCCTGGAATCTCTTGGCAGAGTGGGTTTGTAAACAGGCACCTTTGATCTTCAA
GACGTA CTGTTCCACGGTGGCAGCAGCTTCACAGTTCTGCGGTGCACCAGCTTCGCCCTGGAGAGCTGCG
CCCACCCCGACCCCGCTACTCCCTTGCTGACCTGCTCAAGTACAATTTCTACTTGCCCTTCTTCTTCTT
CGGACCCATCATGACCTTCGACCGCTTCCATGCTCAGGTGAGCCAGGAGCCAGTGAAGCCAGAGGGCGAG
CTATGGCACATTCAGGCCAGGGCGGACTCAGCGCGGCAGCCATCGTGGCTGTTGACGTCTTCTTTCACT
TCTTCTACATCCTCACCATCCCCAGCGATCTCAAGTTCGCCAGCCGCCTCCCAGACAGTGCCTTGTCTGG
CCTGGCGTACTCAAACCTGGTGTATGACTGGGTGAAGGCAGCGGTCTGTTTGGTGTGTCAACACTGTG
GCACGCCTGGACCACCTGGACCCTCCTCAGCCTCCAAGTGCATTACTGCGCTCTATGTCTTCGGGGAAA
CGCACTTCGACCGTGGTATCAATGATTGGCTTTGCAAATATGTGTATGACCACATTGGCGGGGATCACTC
CACCGTATCCCAGAGCTGGCCGCATCTGTGGCCACGTTCTGTCACCACCTCTGTGGCTCGGGCCTTGT
GACATCGTCTACCTGTGGTCTGTTCTTAACTGCTTCGGCCTCAACTTTGAGCTCTGGGTACAGAAGCTGG
CAGAGCGTGGGCCACTGGCACAGATAGAGGCCCGCTATCAGAGCAGATGTCTCGGAGGGTCCGGGCCCT
CTGCGGGGCTGTAAATTTCTGGGCTATCATCATGTACAACCTCGTGAAGCCTGAACAGCCTTGAATTCACA
GAGCTGGTTGCCCGACGCCTGATACTTACAGGCTTCCCCAGACCACACTGGCGGTCTTGTGTTGTACCT
ACTGTGGTGTCCAGCTGGTAAAGGAGCGTGAGAGATCCCTGGCACTGGAGGAGGAGCAGCGGCAGGACAG
AGAGAAGCTGGAGTAGGCATGGGGACAGAGGACCGGGCTCTGCTCAGCTGCTTCTGGGCCCGAGGCCAG
GCCACGGAGCCTGACCAATAGAATAAAAGACTTTTCTACCATAGTAAAAAAAAAAAAAAAAA
```



[View online »](#)

<b>Restriction Sites:</b>	RsrII-NotI
<b>ACCN:</b>	NM_029095
<b>Insert Size:</b>	1512 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">BC024464</a> , <a href="#">AAH24464</a>
<b>RefSeq Size:</b>	1740 bp
<b>RefSeq ORF:</b>	1512 bp
<b>Locus ID:</b>	74770
<b>UniProt ID:</b>	<a href="#">Q9D1G3</a>
<b>Cytogenetics:</b>	9 F4
<b>Gene Summary:</b>	Negatively regulates N-terminal palmitoylation of SHH by HHAT/SKN.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) represents the longer transcript and encodes the functional protein.