

Product datasheet for MC205952

Cyp2g1 (NM_013809) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Cyp2g1 (NM_013809) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Cyp2g1
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC058222

```

GCACCATGATGCTGGGAGGAGCCTTCAGCATCTTTATGGCACTTTGTTTGTCTTGCTGCTCATCCTGAT
TGCTTGGAAGAGGACCAGTAAGGGAGGAAAGCTGCCCCCGGCCACACCAATACCGTTTTTGGGAAAC
TTCTTACAAGTCCGCACTGATGCCACATTTCAAGTCTTCCAGAAGCTCCAGAAAAATATGGCTCTGTGT
TCACTGTGTACTTCGGTCCAAGGCCAGTGGTTGTCTTATGTGGACATGAGGCAGTGAAGGAGGCCCTGGT
GGACCAGGCAGATGACTTCAGTGCCCGTGGTGAATGCCTACACTAGAGAAGAACTTCCAAGGTTATGGT
TTAGCTCTGTCCAATGGAGAACGATGGAAGATTCTTCGGCGTTTTTCCCTGACCGTCCTTCGAACTTTG
GGATGGGAAGCGGAGCATTGAGGAGCGGATCCAGGAGGAGGCTGGCTACCTACTGGAGGAGCTCCATAA
GGTCAAGGAGCTCCCATAGATCCACCTTATACCTGAGCCGCACTGTGTCCAACGTCATCTGTTCCGTG
GTCTTCGGAACGCTTCGACTATCAGGACCAGCGTTTCCAGAGCCTCATGAGGATGATCAATGAGAGCT
TTGTGGAGATGAGCAAGCCCTGGGCACAGTTATATGATATGTAAGTCAAGTCAAGTCAAGTCAAGTCAAG
AAGACACAACCTACCTCTACAACCTTGATAGAAGACCTTAAGGACTTCATTGCATCCAGGGTCAAAATCAAT
GAAGCATCCTTTGACCCCTCAAACCTAGGACTTCATAGATTGCTTTCTCATCAAGATGCACCAGGATA
AAAGTGACCCGCACACAGAATCAACCTCAAGAACTTGTTTCTCACCACCTCAACCTCTTCTTTGCTGG
CACAGAACTGTGAGTTCCACGTTACGGTATGGATTCTTGCTGCTGCTGAAGTACCCGAGGTGGAGGCC
AAGATCCACGAGGAGATTAACCAGGTGATCGGAACACACCGTACCCCAAGAGTGGATGACCGGGCCAAGA
TGCCCTACACTGACGCTGTCATCCATGAGATACAGAGATTGACAGACATTGTGCCCTGGGTGTCCCCCA
TAATGTCACACGGGACACTCATTTCCGAGGCTACCTTCTACCTAAGGGTACAGATGTGTATCCCCTGTTT
GGCTCAGTCCTCAAAGATCCCAAACTTCCGCTACCCAGATGCTTTTTATCCTCAACACTTCTCTGGATG
AACAGGGGCGCTTCAAGAAGAATGATGCCTTTGTGGTCTTTTCTCTGGAAGCGCATCTGCGTTGGGGA
GGCCCTGGCCCGCATGGAACCTTCTCTTACTTCACATCCATTCTTCAGAGGTTCTCCTTACGTTCTCTG
GTGCTCTCTGCTGACATTGACATCGCTCAGAAATCTCAGGCTTTGGCAACATCCCCCAGTTTATGAGC
TCTGCTTCATGGCCCGCTGAGTGCCCACTGCAAGACCAGCAAGAGTCATCTGAGAGCTGAGCTCTGCTG
AAGTTCACACACACACCATCTGCCTTCAATCTCCATACTCCTTGAAGAGAGGTCACTCTCACTGTTG
TTCTGAGACCATGGGTCAAGAGAGGTGTCTGGATATATAAACTAGGGTTCTAGTCACAATCTGGGGGAAA
AGCCTGAGGGTTTCATTAATTTCTGAATTTTCAACCTCAGTCTCAGTGGGCACCAAAATTTACAGATG
CTGTCAGCTTACATAATCTTTGTCTTTGCCCTGATATTGACCAGAATATACCTTGCAATTTCTTCTGTT
TGCTTTTATTTCTCTCCCAATAAAAACTCAAACTAAAAAAAAAAAAAAAAAAAAAAAAAAAAA
  
```


[View online »](#)

Restriction Sites:	RsrII-NotI
ACCN:	NM_013809
Insert Size:	1485 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>BC058222</u> , <u>AAH58222</u>
RefSeq Size:	1886 bp
RefSeq ORF:	1485 bp
Locus ID:	13108
Cytogenetics:	7 A3