

Product datasheet for **MC206003**

Ubl7 (NM_027086) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Ubl7 (NM_027086) Mouse Untagged Clone
Tag: Tag Free
Symbol: Ubl7
Synonyms: 2300004C15Rik
Mammalian Cell Selection: Neomycin
Vector: PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >BC016456
 GGACCAGAGGAGAAAGAGAGAGCGCGAAAGAGAGAGGATGTCTCTCAGACTGGCATCTGGCAGTGAAG
 CTGGCTGACCAGCCACTTGTCCAAAATCTATTCTTCAGTTGCCAGAGACAGAAGCTGGGAGAATATTTCCC
 TAGGGGGCTATAGCATTTCATTTCTGAAGCAGCTTATTGCTGGCAAGCTCCAGGAATCTGTTCCAGACCC
 GGAGCTGATCGATTTGATCTACTGTGGCCGGAAGCTAAAAGATGACCAGACCCCTTACTTCTATGGTATT
 CAGCCTGGATCCACAGTCCATGTTCTCCGAAAATCTTGGCCTGAGCCTGATCAGAAAACCTGAACCTGTGG
 ACAAAGTAGTGCCCTGAGAGAGTTCGGGTCTGCACACTGCCCTGCACAGCAGCTTTCCTACAGAGA
 GGCGGTCTTTAAGATGCTCAGCAATAAAGAATCTCTGGATCAGATCATTGTGGCCACCCAGGCCTCAGC
 AGTGACCCCATTTGCTCTAGGGGTCTCCAGGATAAGGACCTCTTCTGTCTTTGCTGATCCCAACATGC
 TAGATACATTGGTGCCTGCCACCCAGCCCTAGTCAATGCCATCATCCTGGTTCTGCACTCCGTAGCAGG
 CAGCACCCCAATGCCAGGAGCTGACTCCTTCCCGAAGCATGCCCTCCAGCTTTACCGAGATATGCCA
 GGTGGCTTCTGTTGACGGGCTCTCAGATGACGAGGACGACTTTCACCCAAGCACTCGGTCTACGCCCT
 CTAGCAGTACACCTAGCTCCCGCCCGCCTCCCTGGGATACAGTGGAGCTGCAGGGCCCCGCCCATCAC
 CCAGAGTGAGCTGGCCACTGCCCTGGCCCTGGCCAGTACTCCAGAGAGCAGTTCTCATACACCAACTCCT
 GGCACCCAGGGCCATTCTTCAGGCACCTCGCCAATGTCTCTGGTGTCCAGTCAGGGACACCCATCACCA
 ACGATCTTTTAGCCAAGCTCTCCAGCATGCCCTGCAGGCCTCTGGGCAGCCAGCCTTCAGAGCCAGTG
 GCAGCCTCAGTGCAGCAGTTGCGTGACATGGGCATCCAGGATGACGAACTAAGCCTCCGGGCCCTTCAG
 GCCACTGGAGGTGATATCCAAGCAGCCCTGGAACCTCATCTTTGCTGGAGGAGCCCATGAACCTTCTGTA
 CCTGAACACCTAGCAAATTGCAGAGGCTGTTCTGGCAAGAAGCACTTGTGAAGTTGCTTAGATTCCC
 CACCCAACCCACCCCAATATATCATATGCTCCCTTGAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA

Restriction Sites: RsrII-NotI
ACCN: NM_027086
Insert Size: 1143 bp



[View online »](#)

| | |
|-------------------------------|---|
| OTI Disclaimer: | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP). |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | BC016456 , AAH16456 |
| RefSeq Size: | 1340 bp |
| RefSeq ORF: | 1143 bp |
| Locus ID: | 69459 |
| UniProt ID: | Q91W67 |
| Cytogenetics: | 9 |