

Product datasheet for MC206052

Fto (BC057008) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Fto (BC057008) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Fto
Synonyms:	AW743446; mKIAA1752
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC057008

```

GCGGCTTTAGTAGCAGCATGAAGCGCGTCCAGACCCGGAGGAACGAGAGCGGGAAGCTAAGAAAAGTGG
GCTCCTTGAGGAGCTTGAAGACACTTGGCTTCTTACCTGACCCCAAAGATGATGAGTTCTATCAGCAG
TGGCAGCTGAAATACCCTAAACTGGTTTTCCGAGAGGCCGCGCAGCATACCAGAGGAGCTGCATAAGGAGG
TCCCCGAGGCCTTCTCACACTGCATAAGCATGGCTGCTGTTTTCGGGACGTGGTGAGGATCCAAGGCAA
AGATGTGCTCACCCAGTGTCTCGCATCTCATCGGGACCCAGGCTGCACCTACAAGTACTTGAACACC
AGACTCTTACGGTGCCTGGCCGTGAAGGGCTGCACGGTCAAGTACACAGAGGCTGAGATCGCCGCTG
CATGTGACACTTCTAAAGCTCAATGACTACCTCCAGGTGGAGACCATCCAGGCCTTGGAAAGTGGC
TGTCAGAGAGAAGCCAATGAAGACGCTGTGCCACTGTGCATGGCAGAGTTCCCCAGGGCCGGCGTGGGG
CCGTCCTGCGATGATGAAGTGGACCTTAAAGCAGAGCAGCCTACAACGTGACTTTGCTAAACTTATGG
ATCCTCAGAAGATGCCCTACTTGAAGAGGAGCCCTATTTCCGATGGGGAAGATGGCGGTGAGCTGGCA
TACGATGAGAACCTGGTGGACAGGTGAGCCGTGGCAGGTACAGCTATAGCTGCGAAGGCTCTGAGGAT
GAAAGTGAGGACGAGTCCAGCTTGAAGGAGAGATCCTGATACTTGGCATGTTGGTTTTAAGATCTCTT
GGGACATCGAGACACCAGGATTAACAATCCCTCTTACCAGGGAGACTGCTATTTTCATGCTGGATGACCT
CAATGCCACCCACCAGCACTGTGTTTTGGCTGGCTCACAGCCTCGGTTTAGTTCCTCACCAGTGGCA
GAGTGCTCAACAGGCACCTTGGATTATATCTTAGAACGCTGTGAGTTGGCGCTGCAGAATGTCCTCAATG
ACTCAGACGATGGCGACGCTCGTTGAAATCCTTTGATCCTGCAGTTTTGAAACAAGGAGAGGAAATCCA
TAATGAGGTGGAGTTTGGTGGCTGAGGAGTCTGGTTTTCAAGGCAATCGATACAACTTGGCAGGATG
TGGTGGTGTGAGCCCATGACTCACCTGGAGGGCTGTGGAAGAAGATGGAGAGCATGACAAATGGCGTGC
TCCGTGAAGTTAAAAGAGAGGGGCTCCCGTGGAAACAAAGGAGTGAAGTTCTGTCTGCCATCCTGGTCCC
GCTCACCGTGCAGCAACCTGAGGAAGGAGTGGCATGCCAGGTGCCAGTCCCGAGTCCGACTTTTGA
CCAGTACAGCAGAAACCAGACTGCCGGCATAATGGGAGAAGGATGACCTTCCATGCCTCTGCCCTTTG
ACCTCACAGACGTGGTTTCCGAGCTCAGAGGCCAGCTGCTGGAAGCAAGATCCTAGAGGAAGCACAAGTC
TTAGGTGGAGAAGGAAAAGATATCCTTGGCTTTCTTCTCAAACCTTGGACAAGGGTAGTGTAGACTCCT
CTCAACCCAGTGGTTATAGGACCTGGGTCCAATTCAAAACAGCTAGAGAACAAGCCACAGATTCTGA
AGTGTGAAACGTTACTGTGTTTTAGTCTGACTGGTGTAAATTAGAGCACCCCTCCCTCCCCCAAAAT
ATTCTTACTCACCACCCCCACACACACTGGAGCCCTTTGTCTCCCAAGGACAGTAGAGATGGG

```



[View online »](#)

```

TAAGCCTGTGACAAGCCTCACCTGCTCCTCTTTGGCTGGCCACCTCAAAGACCCGGAAGTAGGCCTGG
CCATAGGGTGACAGAAGTAGAGACATCCATCCCAGTTCTCACCCTGGCTTCTAGTGGTTAGACCATCCA
AGATGGCAGCCAGTGTCTCACCGCGGCTGCACAGATATAAATTACAGTCTAGTTGTCATTAACCCATT
CTTCAGTTGCCCTGGCGCATGCTGACTGCTGATTTAGACGACACATACAGATGTTTTATGTGTTCCAC
ATCTGTGAGCGAGCCTTCTTCTCCAGCCTTCTGTAGTGGCTAAAACATGGAGGAAAACAGTAAAAGGC
AGTCACCACTGTACCTTTTCTCTCTGAGAGGGAAGGGCCATGAGCCCTGTCACTGCACACAGGA
AGAGGAAGGAGCTACAAGGAGACAGCCAAACCCTGTCTTTGCGGGGCTTCCCTGATCCAGTGACAGAGAC
CAGCCCTGGGACCACTGTTCTGTGGAGCTAGTGTTCCTCTAGAATTCCCCTCATAGATACTGACAGG
TCGTATGGGGGCTTGAAGCACCCAGTTCCTCAACCCCGTGTGTCAGAGCGTCAGAGAGAGCATGGACTCC
GATTTCAAAGTGCCAGGAGTAAGACCCACCGCAGCACTACTCCATGACCTTTGCATGCTCCTGGGTCTC
AGTTTTCTCATCTGTATGATTTTGGGAGGTACCCAAAATGCCAACTCCTAATGTGAAAATGACGTGATAT
GACAGAGTACAGATCTAGAGTAATTTCTAGCATACCGTGAATATCCGATATGTTCCAGTGTGTTACG
TGAGCTGACCTGGAAGAGGCGTCTTCTCCCCCCCCCCCCCGGTCTGCAGTCTTTGCCGTCTTGT
GTCTCGTGGCCACATCTTTCCCTAGTGACCTGAGGCCTTTTGGTGGCTGAAAGGTTGACAACAGTTAGC
CCATGATTCTCTTGTATTGGCCGAGACCTGCAATTGGAAGGACACAACATAGTTGTGTTTTTGGTGCAT
CCAGGATTCCAGTTAAGGTGCCTTCCTGCCTAGAATCCACCATAGAAGTGCACCTCCACAGCCCAAGG
GCTTCTGGTGCAGGGACCAAGGCTCTCAGACTTTGATTTGTTCCCTTTGCTCTGAGCCAGGACTTC
CTCATCAGTAGAGAGTGAGGTACATTCTCACGGGTAAAATGGAACAACCTAAGTAACGGTATAGGTGGC
TTCATATGTCACAGGCTGGCCCGTCTCCCCATCCCCAGCCATCTGGCTCTGTGGTTTATTAGGCCCGTG
ATATTTCTTTGTTCTATTTAATTGGGTGCCAGAATTAATAATGAGAGATTCCACATAAAATAAAATGT
CTGCATTTCTTTAGAAGCCAGAATACCTGTGAGCCCTGGGTCTCCATAATCGCATGGCCGTGGCTTGG
CAATGCCTGTGACCTTTGTGCTGCTGTGACCCTGTTGCTCCCTGTACCATTTTGTAAACCTGCGTCC
GTCCCAGCAGCACTTTCTAGAGCTGCTTGTGCTGCATCCTTGGTAGAAAGCAGTAGCGTTGAGTGCAGT
CCCTGCTCTCCCATTTGGTCTAAGTTCATGGTCTGCTTGTATTGCTGTGCTGTGATTTGTTTCATACA
CAATAAAGAATTAACGATGAGAAAAA

```

Restriction Sites: Ascl-NotI

ACCN: BC057008

Insert Size: 1509 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC057008](#), [AAH57008](#)

RefSeq Size: 3552 bp

RefSeq ORF: 1509 bp

Locus ID:	26383
Cytogenetics:	8 44.34 cM
Gene Summary:	<p>RNA demethylase that mediates oxidative demethylation of different RNA species, such as mRNAs, tRNAs and snRNAs, and acts as a regulator of fat mass, adipogenesis and energy homeostasis (PubMed:17991826, PubMed:18775698, PubMed:28002401). Specifically demethylates N(6)-methyladenosine (m6A) RNA, the most prevalent internal modification of messenger RNA (mRNA) in higher eukaryotes (PubMed:28002401). M6A demethylation by FTO affects mRNA expression and stability (By similarity). Also able to demethylate m6A in U6 small nuclear RNA (snRNA) (By similarity). Mediates demethylation of N(6),2'-O-dimethyladenosine cap (m6A(m)), by demethylating the N(6)-methyladenosine at the second transcribed position of mRNAs and U6 snRNA (PubMed:28002401). Demethylation of m6A(m) in the 5'-cap by FTO affects mRNA stability by promoting susceptibility to decapping (By similarity). Also acts as a tRNA demethylase by removing N(1)-methyladenine from various tRNAs (By similarity). Has no activity towards 1-methylguanine (By similarity). Has no detectable activity towards double-stranded DNA (By similarity). Also able to repair alkylated DNA and RNA by oxidative demethylation: demethylates single-stranded RNA containing 3-methyluracil, single-stranded DNA containing 3-methylthymine and has low demethylase activity towards single-stranded DNA containing 1-methyladenine or 3-methylcytosine (PubMed:17991826, PubMed:18775698). Ability to repair alkylated DNA and RNA is however unsure in vivo (PubMed:17991826, PubMed:18775698). Involved in the regulation of fat mass, adipogenesis and body weight, thereby contributing to the regulation of body size and body fat accumulation (PubMed:19234441, PubMed:19680540, PubMed:21076408, PubMed:23817550, PubMed:23300482). Involved in the regulation of thermogenesis and the control of adipocyte differentiation into brown or white fat cells (PubMed:19234441, PubMed:19680540). Regulates activity of the dopaminergic midbrain circuitry via its ability to demethylate m6A in mRNAs (PubMed:23817550).[UniProtKB/Swiss-Prot Function]</p>