

Product datasheet for **MC206076**

Scrg1 (BC058532) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Scrg1 (BC058532) Mouse Untagged Clone
Tag: Tag Free
Symbol: Scrg1
Synonyms: AW124307
Mammalian Cell Selection: Neomycin
Vector: PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >BC058532
 AGCCTCCTCGCTTCTCCTTTGCCAGACTGCCCCAGTGAGCTGGAGCATTGAAGAAGAGTCTCCTGCCAA
 TAACACTGAAAAGAAAAGAAAAGGAGCAAGAGCCATGATGAAATCTGTGGTACTTGTATCCTTTGGGCTA
 ACTTTGCTGTTAGAAACACAAGCCATGCCTTCAAGTCGCCTCTCCTGCTACAGAAAAGTTGTAAAGGATC
 GCAATTGTCACAACCTTCCGGAGGGCAGAGCCGACCTGAAGCTGATAGATGCAAATGTCCAGCATATTT
 CTGGGATGGGAAGGGATGCGAGATGATCTGCTACTGCAACTTCAGCGAACTGCTCTGCTGCCAAAAGAT
 GTCTTCTTTGGACCAAAGATCTCCTTTGTGATCCCCTGCAACAACCACTGAGGATCTGCCTTGCACTCTG
 GAGAACATGGTCCTGAAGGCATCTTTGGGCAGCAGAGCCGACCTGAAGCTGATAGATGCAAATGTCCAG
 CATCATTTCTGGGATGGGAAGGGATGCGAGATGATCTGCTACTGCAACTTCAGCGAACTGCTCTGCTGCC
 CAAAAGATGTCTTCTTTGGACCAAAGATCTCCTTTGTGATCCCCTGCAACAACCACTGAGGATCTGCCTT
 GCACTCTGGAGAACATGGTCCTGAAGGCCTTACGTCCCCTAATTTCCCACAACTCTGTCAGTTCAGCG
 CCATTTCTGATATCCATCCAGTATATCCAATCTTCATAGATTCTATAAAGTCTTACTTGCTAGAGTATA
 CTTGGGCTAAAGTGGTAATAAAAGTTGTTTCCATTTGAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA
 AA

Restriction Sites: AscI-NotI

ACCN: BC058532

Insert Size: 541 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online >](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC058532](#)

RefSeq Size: 893 bp

RefSeq ORF: 541 bp

Locus ID: 20284

Cytogenetics: 8 29.8 cM