

Product datasheet for MC206242

2410014A08Rik (BC058245) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: 2410014A08Rik (BC058245) Mouse Untagged Clone
Tag: Tag Free
Symbol: 2410014A08Rik
Synonyms: 2410014A08Rik; AI850530; ESTM19; mKIAA0152; X83352
Mammalian Cell Selection: Neomycin
Vector: PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >BC058245

```
CCGAGGCAGAACGGTCCGCTGCGCTGTTCTCCTGAGTTGGGGCGCGCCGCGGAGGACGAGGGTCCGGCAG
GGGCTGCCCCCGCGGTGGTGGCCCATGCTGAGGCCCGGGAGCCGAGGGAACCGCTGTGGCGTTCT
CCGGCTGCTTCTGTTGCTGCTGCTGCTGGACCGAAGCTCCGGGGTCCGGGGCTCGGCGTGGTCCGGC
GCCGGGGCCGACTGCCGAGAGCGTCATCTGGGCAAGTGAACGCAGGCGCGAGGCGCATGTGGACGTGC
ACGGGATCCACTCCGCAAGGACCCTCTGGAAGCCGGGTGGCCGAGCCTCGGACTATGGCATGAAACT
GCCGATCTGCGCTCCACCCCGGAGGACCAGATCCTATATCAGACAGAGCGGTACAATGAAGAGACCTTT
GGCTACGAAGTACCCGTCAAGGAGGAGGGCGACTATGTGCTGGTGTGAAGTTTGTGAGGTCTACTTTG
CACAGTCTCAGCAGAAGGTGTTGATGTCCGTTGAATGGCCATGTGGTGGTGAAGGACTTGGACATCTT
TGACCGTGTCCGGCACAGCACGGCCCATGATGAAATCATCCCGATGAGCATCAGGAAGGGGAAGCTGAGT
GTGCGGGGGGAGGTGTCCACCTTACAGGGAAGCTATACATCGAGTTTGTCAAGGGGTACTATGACAATC
CCAAAGTCTGTGCACTCTACATCCTGGCTGGGACAGTAGACGATGTACCAAAGCTGCAGCCTCATCCAGG
ATTGAAAAGAAAAGAGAGGAGGAGGAGGAGGAAATGATGAGGGGTCTAATCTCAAAAGACAGACC
AACAAAGACAGGGTGCAGTCAGGCCCGCACGCCCTAACCCCTATGCCTCGGACAACAGCAGCCTCATGT
TCCCTATCCTGGTGGCCTTCGGAGTGTTCAATCCAACCCCTCTTCTGCCTCTGCCGCTTGTGAGAGGACAA
ATGGCCATCCTCAGCAGGGTGGAGGGAGTGTGCGAGAAACCAGACACTTTGGTTTCTGTAATCTTACAG
TGTTTAGGAAAAGAGAAAACAAATCAAACAGCGCACACACACACACACACACACACACACACACACACAC
ACACAAACACAGACACACATGTGCTCACACACACACGGAAGAGATAAGAAAAGCAGAGAGAGTGGGGTC
ACCCACCAACACAGGTTGCAGAGTGCTCAGTCTGCTGCTGCTAGCTGGCCCCACCCAGCTCTCTTTCT
TCTGGGGTACTTAGGGAGAGGATCCTTCTCCTCCGAGGATCCTCATTATATACCTAGTAGGAAGGACTCT
GAACTCAGAGCAGTCACTTCTTGGTTTTTTAGGTTAGACATTAACAAACAGCAGGAAAATGTCACA
TAACCTGGAAAGACCGCACTGGAGTTGGGTGTGTTCTGAGTTTCTGTCCAGATGTGATTTAGCAGTCAT
GGTATTCCAAGTGTGTGTGGAGGTAGATAAGTGGGACGGTAAAGACAAAGGCGTGGAAGAGAAAGTCA
GCGCTTCTTCTGGGAAGCCTGTCGTCTGGGTTAATCTGTTGATGTATCTTATTGTCATCTCTCATCC
TGATGTGTGCCAGTTAGAGGACTGTCTGGTATCCAAGACGACTAGTGTATGTCAGGTCCTCAAATCACCA
ACAGCTTGTTACTACAACAAATTGTTTTGATTTTAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA
AAAAAAAAAAAAAAAAAAAAAAAAAAAA
```



[View online »](#)

Restriction Sites:	RsrII-NotI
ACCN:	BC058245
Insert Size:	876 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC058245 , AAH58245
RefSeq Size:	1772 bp
RefSeq ORF:	876 bp
Locus ID:	109154
Cytogenetics:	5 55.99 cM
Gene Summary:	Carbohydrate-binding protein with a strong ligand preference for Glc2-N-glycan. May play a role in the early steps of protein N-glycosylation (By similarity).[UniProtKB/Swiss-Prot Function]