

Product datasheet for **MC206274**

Rpl27 (BC082284) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Rpl27 (BC082284) Mouse Untagged Clone
Tag: Tag Free
Symbol: Rpl27
Mammalian Cell Selection: Neomycin
Vector: PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >BC082284
GCTGGGTGCCGGCTGCCTGCTGTCGAGATGGGCAAGTTCATGAAACCCGGGAAAGTGGTGTGGTCCTGG
CCGGACGCTACTCCGGACGCAAAGCCGTCATCGTGAAGAACATTGACGATGGCACCTCGGACCGCCCTTA
CAGCCATGCCCTGGTGGCTGGAATTGACCGCTATCCCCGGAAAGTGACAGCTGCCATGGGCAAGAAGAAG
ATCGCCAAGCGATCCAAGATCAAGTCCTTTGTGAAAGTTATAACTACAACCACCTCATGCCACAAGGT
ACTCTGTGGATATCCCCTTGGACAAGACTGTTGTCAACAAGGATGTGTTCAAGGACCCAGCTTTGAAGCG
CAAGGCCAGGGCGGGAGGCCAAGGTCAAGTTTGAGGAGCGATAACAAGACAGGGAAGAACAATGGTTTTTC
CAGAAGCTTCGCTTTTAGGTATATTTTTGTTTTGGTCATTAATAAATTAATAAAAAAAAAAAAAAAAAAAAAA A

Restriction Sites: EcoRI-NotI

ACCN: BC082284

Insert Size: 411 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC082284](#), [AAH82284](#)

RefSeq Size: 491 bp

RefSeq ORF: 411 bp

Locus ID: 19942

Cytogenetics: 11 D

Gene Summary: Component of the large ribosomal subunit (By similarity). Required for proper rRNA processing and maturation of 28S and 5.8S rRNAs (By similarity).[UniProtKB/Swiss-Prot Function]