

## Product datasheet for **MC206383**

### Eif1b (BC033505) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Eif1b (BC033505) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Eif1b  
**Synonyms:** 1500010M16Rik; 3110001N14Rik; Gc20  
**Mammalian Cell Selection:** Neomycin  
**Vector:** PCMV6-Kan/Neo (PCMV6KN)  
**E. coli Selection:** Kanamycin (25 ug/mL)

**Fully Sequenced ORF:** >BC033505  
 GAAGCCGCAGCGCCCTCCCTGCTCGCGCCCCACCGCCGCTCCGCGTCTTCTGCGCCCCGGGCTTCC  
 GCCGCTTCATTGATTTCCGTTTTCTCGCGCTGCAGCCTCCTGACACGGTGATCCGGGCGGGCCCCGGGG  
 AATTCTATCCTCTCACCGCCTCACACTAGTGTGCGCATGTCCACTATCCAGAACCTCCAATCTTTGACC  
 CCTTTGCTGATGCAACTAAGGGCGACGACTTACTCCCGGAGGGACTGAGGACTACATTCATATAAGGAT  
 CCAGCAGCGCAACGGCAGGAAGACGCTGACCACTGTGCAGGGCATTGCGGACGATTATGACAAAAAGAAA  
 CTTGTGAAAGCTTTCAAAAAGAAATTTGCCTGCAATGGAAGTGTGATTGAACACCCTGAGTACGGAGAGG  
 TTATTCAGCTTCAAGGGGACCAAGGAAGAATCTGCCAGTTTCTTCTGGAGGTTGGCATCGTCAAGGA  
 GGAGCAGCTGAAGGTTACGGGTTCTAAGATGAACCAAGCGTGGCGAGTTTCTTAAATGGTTTTGTCT  
 CTAAGTCAAGTTGCTGCTCGGAGATGATTCTTTACAGTAAACGACAGACTTTGCGTTTTATTTAATCA  
 TTCAGACTTCCACTCACACTGCATGGCTACAGAAAACACGGGGTATGTAGGCTCCTGAGTCATAAGGAA  
 ATCGCGGTGAGATGGGAACGAAGCCGAGTCCGTCTAATGTTTCCAATGGAAAAGGTTTTGTTGTGA  
 GTGTTTGTGTTAGTTTTATTTCACTTACTGTTAAATGTTTTGTTGTTGATTAAACCATGTATGTT  
 GCGGCTTAAACAATAA

**Restriction Sites:** RsrII-NotI

**ACCN:** BC033505

**Insert Size:** 342 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [BC033505](#), [AAH33505](#)

**RefSeq Size:** 910 bp

**RefSeq ORF:** 342 bp

**Locus ID:** 68969

**Cytogenetics:** 9 F4

**Gene Summary:** Probably involved in translation.[UniProtKB/Swiss-Prot Function]