

Product datasheet for **MC206406**

Olfml1 (BC047207) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Olfml1 (BC047207) Mouse Untagged Clone
Tag: Tag Free
Symbol: Olfml1
Synonyms: MGC56882, MVAL564, 6720478C22, ONT2, mONT2
Mammalian Cell Selection: Neomycin
Vector: PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >BC047207

```

GGAAGCTGTAAGTAGCCATTAGCCTGGAGAGGGTCTGTGCATGGCCAAGCTGAGGGCAGCAGACACTGGA
CCCCGGAGCTCTCTGCAACCAAGTGAACACGTAAAGCCTGCAGAGACTGGCTTGCCAACAGCTGGACTT
GATCACCAGCTTGAAACTGAGATCATGCACCGGTAGAATAACAAGCGGACTTTGTTCTCTCTGGCTATG
CAAAGCGTTGTTTTCCAAGGAATTGCCACAGCAGCAGACAGAAAACCAGAGATCACCACCCGGACCCTG
AGTCGTTTTAGAGAACCAGAGTCTCGGGAAGTGTGCAGGATGGTGGTGGCCCTCAAGAAGCTTGTCT
TCCCTGGTTCTCTTCTTGCAGCATTCTGCCCCACCACAACATGCCCAGGACCCAGCCATGGTACACT
ACATCTACCAGCGCTTCCAAGTCTCGAGCAAGGGCTGGAAAAATGTGCCCAAACAAGGGCATATAT
TCAAGACTTCCAGGAATTCTCAAAAACATATCCATCATGCTGGGGAGGTGTGAGCCACACAAGTGAA
TATAAGAGTGCAGTGGAAAACCTTGCCCTAAGAGTGGAGCGTGCCAGCAGGAGATCGACTACCTGCAGT
ACCTGCGAGAAGCTGACTTCTGCATTGAATCTGAGGAGAAGACTGGCTGAAAAGTTGCTTCAAGAAGA
AGCAGAAGAAAAGATCCGAACCCTGTTGAACAACAAGCTGTGACAACATGCTGATGGCTATAAAGTCT
CTGAAAATAGTGAAGAAGACCGTGGACCCAGATGGCTCTTGGATGAAAGATGCTGGCAGTAACTCTGCAA
AAGTGTATCTATTAGCTGGATCCAGAAAACAACAGTGGGAAATTTGCAAACCTGCGGGCATTCTATGGA
AGATAGCATCAAGCCTGGTCCCCGTAAGTTAATCTTACCCTTCTGGCAGGGATCAGGGCAAGTGGTC
TACCAAAGCTTTCTATTTTTTTCACAATCAAGGAACCTTAAATGAGATAATTAATAATACTGCAGAAG
AAGACTGTGGAAGATCGAATGCTTCTCCAGGAGGGGCAGGCCGAGCACCATCTATCAACTCTCTCT
CTACGTACATTGACCTTGTGTGGATGAACACGGACTTGGGCCATTCAGTACAGGATACCCGTGTAGA
CCATTTGGTTCTCACAAAATTGAGCCTGACACTTTGGGAGTAGAACATTCATGGGATACCCGTGTAGA
AGCCAGGATGCGGAAGCTTCACTTCTACTGTGTGGGTCTCTATGTTGTCTACAGTTCGGTGGCCAGG
GCTCTCATCGAATCACCTGTGTCTATGATCCACTGGGCACTGTGAGGAGGAGCATCTGCCCAATTTGTT
CTTCCCTAGGCGGCCAAGAAGTCACTCCATGATCCATTACAACCCAGAGATAAGCAGCTTTATGCCTGG
AATGAAGGCAATCAGATCATTTACAACTCCAGACAAAAGAAAAGCTGCCTCTGGAGTAATGGCTGGCAT
TGGAGGGAGAGCAGTGTGACTTGGGCAGTGGTTTTCCGGGACACTGAGGTACAGCCCTTTATAGTGA
GTATCTGGCCCTCTCGATGCACATAGGAATAGGTTCTAAAAGTTGAAATACATATTTCTCTTTCCCAA
TGTCAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA

```



[View online »](#)

Restriction Sites:	RsrII-NotI
ACCN:	BC047207
Insert Size:	702 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>BC047207</u> , <u>AAH47207</u>
RefSeq Size:	1748 bp
RefSeq ORF:	702 bp
Locus ID:	244198
Cytogenetics:	7 E3