

Product datasheet for **MC206492**

BC008163 (BC008163) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	BC008163 (BC008163) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	BC008163
Synonyms:	MGC11672
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for BC008163, the custom clone sequence may differ by one or more nucleotides

```
ATGGCGGCGCTGTACGCCTGCACCAAGTGCCACCAGCGCTTCCCCTTCGAGGCGCTGTCTCAGGGGCAGC
AGCTGTGCAAGGAATGCCGATTGCTCACCTGTTGTGAAGTGTACCTACTGTAGGACCGAGTACCAGCA
GGAGAGTAAAACCAATACAATTTGCAAGAAGTGTGCTCAGAATGTACAGTTGTATGGAACGCCTAAGCCT
TGTCAGTATTGCAACATAATTGCTGCATTTATTGGGAATAAATGCCAGCGCTGCACAAATTCAGAGAAGA
AGTATGGACCCCATATTCCTGTGAACAGTGAAGCAGCAGTGTGCATTTGACAGGAAGGATGATCGAAA
GAAGGTGGATGGGAAGTTGCTATGCTGGCTGTGCACGCTGTCGTACAAGCGTGCCTTCAGAAGACCAA
GAGCAGAGGAAGCACCTGAGCAGCTTCTCGAGGCAGCCATCAGGAGAAGGAGCAGTACAGCCGCCTGA
GTGGAGGCAGCCATTACAACAGCCAGAAAACACTTTCTACATCTTCAATTCAAAATGAAATCCCAAAGAA
AAAATCCAAGTTTGAGTCAATCACAACCTAACGGAGACAGCTTTTCTCCAGACCTGGCTCTGGATTCCCA
GGCACTGACCACTTTGTGTCATCATTGCTCAGCTGAAGGAAGAGGTGGCCACTCTGAAGAAGATGCTGCATC
AGAAGGATCAAATGATCCTAGAGAAAGAGAAGAAGATCACGGAGCTGAAGGCTGATTTTCAGTACCAGGA
GTCTCAGACTCGGGCCAAAATGAACCAATGGAGAAAACCCACAAGGAGGTCACAGAGCAACTGCAGGCC
AAAAACCGAGAGCTTCTGAAACAGGCAGCTGCCTTGTCCAAGAGCAAGAAGTCTGAGAAGTCAGGCACCA
TAACGTCTCCGTGA
```

Restriction Sites:	EcoRI-NotI
ACCN:	BC008163
Insert Size:	924 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC008163 , AAH08163
RefSeq Size:	1615 bp
RefSeq ORF:	924 bp
Locus ID:	230789
Cytogenetics:	4 D2.3