

## Product datasheet for **MC206537**

### Mrpl35 (BC028750) Mouse Untagged Clone

#### Product data:

|                                  |                         |
|----------------------------------|-------------------------|
| <b>Product Type:</b>             | Expression Plasmids     |
| <b>Tag:</b>                      | Tag Free                |
| <b>Symbol:</b>                   | Mrpl35                  |
| <b>Synonyms:</b>                 | 1110066C01Rik           |
| <b>Mammalian Cell Selection:</b> | Neomycin                |
| <b>Vector:</b>                   | PCMV6-Kan/Neo (PCMV6KN) |
| <b>E. coli Selection:</b>        | Kanamycin (25 ug/mL)    |

**Fully Sequenced ORF:** >BC028750  
 CCGTAGGCAGGCTGCAAGATGGCGGCTTCCATCTTCACGGGCGCCGTGAGAGCTGCCTCAGGAATATTTT  
 GTCCCTTGAATGTCTTGGCATCTTCAACCTATCGGAAGTGCAGCAAGTGCCTGTCTTAATCTTCACT  
 GTGCACCATTCATTTAGACATATTCAGACATCAGTTGTTTCTCTGCTCCAAGACTTGTCACATCTGTT  
 GGACACCTGGCATATGGGCATACCACAACAGTCCTTAATAGAGTGGCCACATTGGTCCCAGTGCCTGA  
 AGCCACCAAGTCAGAGCTTAACATACTGCAGCACACGAAAGGCAAGAGAAAGACGGTAAAGTCGGTCGT  
 GCATCGCTTCTCCGACTTCATTCTGGCCTTTGGCTAAGGAGAAAGGCTGGCTATAAGAAAAAATTATGG  
 AAAAAGTCAACTGCAAGAAAGAAGCGCTTGAGGGAATTTGTGTTCTGCAGCAAGACCCAGAGTAACTCC  
 TGGATAAAATGACGACATCTTCTGGAAGAGGCGCAACTGGTACGCTGGCGATCCTTACCAGATGTACCA  
 CGATCGAACAACCTACGAGTATAGTTTACAGCTATCATAGTTCCAGCTATGGAATGTGGGCTTTGTG  
 TACATAATGCCCTTGCAAGTGAGTATGCGTTAACTGTAAGTGTACAGCTGATCGCAAGCTTAGAT  
 TACCACCATTTGTGATTAAACATTTCCAGAGTTGGTCATAGAACTAAAAGTTAAATTTATTGTGATTAT  
 TTCAAAAACTGGTAATGGCTTTCTTATTTCTGGAATCTTGGATGAGGAAAAACACGAAAGACGCTCA  
 ACAGTGGGCTTTGTGGAATTGAGATGGTATCCAGGTTTTTGCACATGCTAACAGACATCTTACATG  
 CCTATTTCCCAATATGACACATGCTGTGTTAAAAAATAAATACCTAGGAAGTACAGAAAAAATATATATAT  
 AAAAAA

**Restriction Sites:** EcoRI-NotI

**ACCN:** BC028750

**Insert Size:** 567 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



|                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).                                                                                                                                                                                                                                                                                                                                    |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol> |
| <b>Note:</b>                  | Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>RefSeq:</b>                | <u><a href="#">BC028750</a></u> , <u><a href="#">AAH28750</a></u>                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>RefSeq Size:</b>           | 987 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>RefSeq ORF:</b>            | 567 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Locus ID:</b>              | 66223                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Cytogenetics:</b>          | 6 C1                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |