

## Product datasheet for **MC206566**

### Nkd1 (BC034838) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Nkd1 (BC034838) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Nkd1  
**Synonyms:** 2810434J10Rik; 9030215G15Rik; Nkd  
**Mammalian Cell Selection:** Neomycin  
**Vector:** PCMV6-Kan/Neo (PCMV6KN)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Fully Sequenced ORF:** >BC034838

```

CGACGGCGGCAGGAGCGCGTCCCGGCCCGCTAAGGCTGCGCTCGGCGCGGACTGTGAGGAGGAGGC
GAGCGAGGCTGGCGCGGGTTCGGCGGCCGACGCATGGCTTAGGACGCTCCGCCCGCGCCCCAGCATG
GGGAAACTTCACTCGAAGCCGGCCGCGTGTGCAAGCGCAGGGAGAGCCGGAAGGTGACAGCTTTGCTG
TAAGCGCTGCTTGGGCAAGGAAAGGCATCGAGGAGTGGATCGGGAGGCAGCGCTGTCCAGGCAGCGTCTC
AGGACCCCGTCAGCTGAGATTGGCAGGCACTGTTGGTCGAGGCACTCGGAACTCGTGGGTGACACTTCT
AGAGAGGCTCTCGGTGAGGAGGACGAGGACGACTTCCCCCTAGAAGTGGCCCTGCCCGCTGAGAAGATTG
ACAGCCTAGGTAGTGGAGATGAGAAGAGAATGGAGAGACTGAGCGAACCTGGCCAGGCTCCAAGAAGCA
GCTCAAGTTTGAAGAGCTGCAGTGTGATGTCTCTGTGGAGGAGGACGCCGGAAGAGTGGACTTTCACT
CTATATGACTTCGACAACAATGGCAAAGTGACCCGTGAGGACATTACCAGCTTGTGCATACCATCTATG
AAGTGGTTGACTCCTCTGTGAACATTCCCCACATCAAGCAAGCACTGCGGGTGAAGCTCACCGTGGC
TCCTGACGGGAGCCAGAGTAAGAGGAGCGTCTTTTCAACCATACCGATCTGCAGAGCACAAAGCCCGCA
GCAGACACCAAACCCGCTGAGGAGCTGCGTGGCTGGGAGAAGAAGCAGCGAGCCCACTCAGGTTCCAGG
GTGACAGCCACCTGGAGCAGCCAGACTGTACACCATTGCGTGGATGAGAACATTGAGAGGAGAAAACCA
CTACCTAGACCTGGCGGGATAGAGAACTACACGTCTCAGTTTGGACCGGATCCCTTCGGTGGCCAG
AAGTCAGAGCTGCCCCCTCGAATCTCCAACCCACTCGCTCTCGCTCCCACGAGCCAGAAGCTGCCACG
TCCCACACCGGAGGCCCAAGGTGTGGACCCAGGCTCCTTCCACCTCCTTGACACCCATTTGCCAAGGC
ATCAGAGCTCCAGCAACGGCTCCGGGGCACTCAGATGGGAGCAAGCACTTTGTGAGTCCCCAAGGCC
CAGGGCAAGAACATGGGTATGGGCCACGGGGCCAGAGGTGCAAGAAGCAAGCCCTCACTGGTACCCACCA
CCATACTGTCTCCCCCTCTGCCATCTGGCCACCAAGCCAGCCCTTCTCCCCACCTGGCACCCCTGGG
GCACAAGAAACACAAGCATCGAGCCAAGGAGAGCCAGGGAGCTGCCGGGCCTGCAGGGCCCCCTGGCT
GCAGGAGGCTCCACCGTCATGGGGCGGAGCAGGTGAGGGAGCTGCCTGCCGTGGTGGTGTACGAGAGCC
AGGCTGGGCAGGCCGTCCAGAGACACGAACACCATCACCACCAGAACATCACCACCATTACCACCATT
CTATCAGCCCTAGACCCAGCAGGCTGCCACGGGAAGGACCCAGCCACACCCCTAAGGCATTATTATTCT
ATTAATTATTGTTATTATGGCAATTATTGTTATTAATAATTATTGTTACTCCACTAATATTTAGCCAGCC
TTCATGTAGAAGACACATGGAAACACAGAAGTAACTTTTATTTATAAAAAAAAAAAAAAAAAAAAAAAAAA
AAAAAAAAAAAAAAAA

```



[View online »](#)

<b>Restriction Sites:</b>	EcoRI-NotI
<b>ACCN:</b>	BC034838
<b>Insert Size:</b>	1416 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">BC034838</a> , <a href="#">AAH34838</a>
<b>RefSeq Size:</b>	1765 bp
<b>RefSeq ORF:</b>	1416 bp
<b>Locus ID:</b>	93960
<b>Cytogenetics:</b>	8 C3
<b>Gene Summary:</b>	Cell autonomous antagonist of the canonical Wnt signaling pathway. May activate a second Wnt signaling pathway that controls planar cell polarity. Required for spermatogenesis. [UniProtKB/Swiss-Prot Function]