

Product datasheet for **MC206619**

Ppic (BC025861) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Ppic (BC025861) Mouse Untagged Clone
Tag: Tag Free
Symbol: Ppic
Synonyms: CyP-20c
Mammalian Cell Selection: Neomycin
Vector: PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >BC025861
 CGTCCGGCCACAGTCCCAGGAGCCTAAGGTTGCGCCCGTGCCATCGTTGCCACTTCCACCGCGCACCATGA
 GCCCGGGTCCCCGTCTGCTGCTGCCCGGGTGCTCTGCCTGGGGCTTGGAGCCCTGGTGTCTTCTTCGGG
 GTCCCTCAGGCGTCCGAAAACGAGGTCCTCGGTGACGGACAAGGCTTCTTTGATGTGAGGATCGGAGAC
 AAAGATGTGGGTAGAATTGTGATTGGCCCTTTTGAAAACGTTGTACCCAAGACGGTGAAAACTTCGTGG
 CTCTGGCAACAGGAGAGAAAAGGCTACGGGTACAAGGGCAGCATCTTCCACCGTGCATCAAGGACTTCAT
 GATTCAAGGAGGGGACTTCACAGCTAGAGATGGCACTGGAGGGATGAGCATCTATGGTGAGACATTTCCG
 GATGAGAACTTCAAAGTGAAGCACTATGGCATTGGCTGGGTGAGCATGGCAATGCTGGACCAGACACCA
 ACGGATCCAGTCTTTCATCACCTTGACCAAGCCACCTGGTTGGATGGCAAACAGTGTATTTGGAAA
 AGTCCTAGATGGGATGACTGTGGTACATTCCATTGAACCTCAAGCAACAGATGGCCATGACCGCCACTC
 ACCGACTGCACCATCGTCAACAGTGGCAAGATAGATGTGAAAACACCCTTGTGGTTGAGGTCCCTGATT
 GGTGACAGAAATGGCGGAAGACAGGAAACACACTGGCTGGCCCATGTGCACATCTGTGATCCACCAGCAG
 TCTTTTAGTTACTTAAAAGTTGTAGTTAAAAAAAAAAAAAAAA

Restriction Sites: EcoRI-NotI

ACCN: BC025861

Insert Size: 639 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC025861](#), [AAH25861](#)

RefSeq Size: 812 bp

RefSeq ORF: 639 bp

Locus ID: 19038

Cytogenetics: 18 D1

Gene Summary: PPlase that catalyzes the cis-trans isomerization of proline imidic peptide bonds in oligopeptides and may therefore assist protein folding.[UniProtKB/Swiss-Prot Function]