

## Product datasheet for **MC206707**

### Upk1b (BC046604) Mouse Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Upk1b (BC046604) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Upk1b
Synonyms:	UPIb, Tspan20
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC046604

```

CGGAGAAGAAGGCAGAGGAGACCAGGCAGCCGGTCTTTTAGAAAATCCTGAAGATGGCCAAAGACGATTC
CACTGTTTCGTTGCTTCCAGGGCCTGTTGATTTTTGGACACGTAATTGTTGGTATGTGTGGCATCGCCCTG
ACAGCAGAGTGTATCTTCTTTGTATCTGACCAACACAGTCTTTACCCACTTCTCGAAGCCACCAACAATG
ATGATATCTTCGGGGCAGCTTGGATAGGCATGTTCTGGGCATCTGCCTCTTCTGCTTGTCCGTTCTAGC
CATAGTAGGAATTATGAAGTCCAACAGGAAAATCCTCTTGGCGTACTTCATCATGATGTTTATAGTGTAC
GGTTTTGAAGTGGCATCTTGCATCACAGCAGCAACACAACGCGACTTTTTTACAACCAACCTCTTCTGA
AGCAGATGCTGATGAGGTATCAAAAACAGTCCCCAACCAATGACGACGAATGGAAGAACAACGGGGT
CACCAAGACCTGGGATAGGCTCATGCTGCAGGACCACTGCTGTGGTGTAAACGGTCCGTCAGACTGGCAG
AAATACACCTCTGCCTTCCGGGTGGAGAATAACGATGCTGACTACCCCTGGCCTCGGCAGTCTGCGTCA
TGGACAAGCTTAAAGAGCCTCTCAACCTGGACGCTTGCAAACCTGGAGTGCCTGTTATTACCACAGTCA
GGGCTGCTATGAACTGATCTCTGGACCAATGGACCGGCATGCCTGGGGAGTTGCCTGGTTTGGATTTGCC
ATCCTCTGCTGGACTTTTTGGGTCTCCTGGGCACCATGTTCTACTGGAGCAGAATTGAATATTAAGCAC
AGACAGTGTCTCCACCATACCTCCCTCCCTCCCTGGGTGACTCTGAATTTGGTGTGGAACCTGCTGTTT
CGTGACACGCCCTCTGTGTGCTCTGCCCGTGACTCGTGCTTGTCTCTCAAATGCCAGTTCGGGTGGA
ATGGGCTTCCTTAGGCTTTCAAGATGCTGCAATCCTTAGCAGAAATATTTCTCATCCTTGCCCTAAGATT
CCGAAATGATTTATACTATAAGACAGAAGGTGACTATATATGACTATATAATATAGTACTTAC
GTGACTATATAAGACAGAAAGTTTGGAGAATGGAGACGACTCTCATTCTGTCTATGTTAATGTGGAG
GCACAATGATATTCAAAGAACCTGCTTACCATTATTATAAGCTCTGTCTGTATTGAGCAGCAAATTC
AAGAAACAGAGAGCCCATCTCTATTTTCTCAACCTTGAAAGACTTTTGTGTTGCACTGTGTAGGCTG
GGGGAGGAGGGTTGCACTGTGAAAGAGGAGGAAATGGCAGTCATGCTTTCTTCTGCCTCTGGTTCTGCTG
TGACTTCTGTCTGCAGTGTCTCACACTTTTGGATGCTTTTAGTCTCCAGTACAACCTTGATTTTCAG
ATTTCTCTCATGCCCTGTGGAGTGTGTGATTGAGTGTATGGAACTTTGGGGGTGTCGGCCGTGCC
GATTTACATCAGCTCAGTCTTTGCCCTCAGAGTTTCAACCGGAAACACTTCCCAAACCGCGGGATCCAT
GCCAGCTATTCAAACAAGTGATTTATAGGGTCTTTTGTGCTGTCAAAATAGTTTGTCTCACATTAACGA
TCCTGAAGGAAGTCTGACAAGTACCATTGTTTGCCTATCCATAAATAAAGGAGCCCATAAAAATAAAA
AAAAAAAAAAAAAAAAAAAAAAAAA

```



[View online >](#)

<b>Restriction Sites:</b>	EcoRI-NotI
<b>ACCN:</b>	BC046604
<b>Insert Size:</b>	783 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">BC046604</a></u> , <u><a href="#">AAH46604</a></u>
<b>RefSeq Size:</b>	1773 bp
<b>RefSeq ORF:</b>	783 bp
<b>Locus ID:</b>	22268
<b>Cytogenetics:</b>	16 B4
<b>Gene Summary:</b>	Component of the asymmetric unit membrane (AUM); a highly specialized biomembrane elaborated by terminally differentiated urothelial cells. May play an important role in normal bladder epithelial physiology, possibly in regulating membrane permeability of superficial umbrella cells or in stabilizing the apical membrane through AUM/cytoskeletal interactions (By similarity).[UniProtKB/Swiss-Prot Function]