

Product datasheet for MC206742

S100a11 (BC021916) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	S100a11 (BC021916) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	S100a11
Synonyms:	EMAPI, cal, S100c, S100a14
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC021916 CAGCCACAGTCGCAGCCGACCTCGCTCCTCAACTTGAAGCAAAATGCCTACAGAGACTGAGAGATGCA TTGAGTCCCTGATTGCTGTTTTCCAAAAGTACAGCGGGAAGGATGAAACAACACTCAACTCTCCAAAAC TGAATTCCTTTCTTCATGAACACAGAGCTGGCTGCCTTCACAAAGAACCAGAAGGATCCTGGTGTCTT GACCGCATGATGAAGAAGCTGGACCTCAACTGTGACGGGAGCTAGATTTCCAAGAGTTTCTCAACCTCA TTGGTGGCTTAGCTATAGCGTGCCATGATTCTTTCATCCAAAATTCCCAGAAGCGAATCTAATCCTCTTA GTTTCCATCCAACCACCAAGTCATCACCTCCCCAATCCACACCTGCACTGAGCCCAGCACACCT ACCACAACATGTATCCCACGCCTGCTGGCGAAAAATAAACAAATGCCATTTTTTTTTAAATGTAAAAAAA AA AAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA
Restriction Sites:	RsrII-NotI
ACCN:	BC021916
Insert Size:	297 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC021916](#), [AAH21916](#)

RefSeq Size: 585 bp

RefSeq ORF: 297 bp

Locus ID: 20195

Cytogenetics: 3 F2.1

Gene Summary: Facilitates the differentiation and the cornification of keratinocytes.[UniProtKB/Swiss-Prot Function]