

Product datasheet for MC206743

Clec2h (BC021766) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Clec2h (BC021766) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Clec2h
Synonyms:	clr-f; Clrf
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC021766

```

AGTTGAGTCATAGTTCTTCTCTCAGCTCACAGATGAATGCTGCAAAGGTTGAAATGAGTTCATGGGCA
TGCTGCAGAGGGCAGACCTCACTGCTGCAGACTGCCTGCAGGAGGAGAAATGGGTAAAAAATCCAAGG
AAAATGTTTCAGAATCATCTCCACTGTGTCTCCTGTTAACTTTACTGCTGCTATGGAGTGATCATGGTC
CTCACTGTAGCTGTAATTGCTCTTTCACTTGTCTTGTGCTGAGAGAAATAAAATCCAGCCATGGAAGACC
GTGAACCTTGCTACACCGCATGCCAAGTGGCTGGATTGGATTGGAAGTAAATGCTTTTATTTTCTGA
AGACATGGGAAACTGGACCTTCAGTCAGTCCTCTGCGTAGCATCAAAATCCCATCTAGCACTGTTTCAC
AGCCTGGAGGAGCTGAATTTTCTGAAGAGATAAAGGGGACTTCTGACCACTGGATTGGCTGCACAGAG
CATCAACACAGCACCTTGGATATGGACAGACAACACTGAATATAGCAACTTGGTTCTCACGAGAGGAGG
AGGAGAATGTGGCTTCTGAGTGATAATGGGATCAGCAGTGGCAGGAGCTACACACATAGGAAGTGGATT
TGTAAGCAAGTTTGTGAGTAGTTGTAATCCAGAGTAGGAAGTGTACCAAGGCATGTGTGAGAATCAATCA
GATGCTCTATGACCTGTTATCTACAGTTGGTGCTCTGCTCTGATCAACCAAAACCATTATCAAAAACCA
AAACCATTATCAAAATCCAGCAATAATCACATTGACAGCCAAGATGAAAATTAGTTTAGAGGAGTGTGG
GCACCGGACAACCTGAGTAGGAGACAACCTGAGTCCCGCCAGAGTACTGGTGCTCTGGGCGGCCATACTC
AGGAGTACTGTGCAAACTTTCAACGTGGCCCCGGGTGGGTGTCTGGCTGTGGGAAGCCAGAGAACCCCA
GCCCCGCGGAGTTGGGGTGTGTGGACAAGTGGCATTCCGAGGTCTCTGAGCCTCATGGATGCCAGAGATA
ACAGGTTCTCCGTCTCGAGCATCCCAGACACCGTTGGACACTGTGAAAGAGCTAAGGACAAAGTGGGAGC
CCAGGGGGCAGCTTGGACAATTTTCTGATCTGAAAAGAATCTAGTAAAGAGTAAATTCAAATGGCAT
ACACACGTATGGAATTCTCAAAAAGATGAGACTCTCGCCTCGCCCCCTGCTGTCTCCTCCCTTTGGTCC
TGGGGATGGCCAAGCACTATTTATTGATATTGATCTTACTTCTGTGCTAAGCACATTATTTTGAGAAAT
GTTTTGCAGTTCTTATCAAGTTTGTGTATCTGTTTTTACTGAAGTCTTTGGTCCATGTCGAACCTAAG
TTTGGTACAGAGTGATAAGTATGATATATTTGCATTTTCTAATAAAGACAGTTTGAAGATCAAAAAA
AAAAAAAAAAAAAAAAAAAA
  
```

Restriction Sites:	RsrII-NotI
ACCN:	BC021766


[View online »](#)

Insert Size:	657 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	BC021766 , AAH21766
RefSeq Size:	1490 bp
RefSeq ORF:	657 bp
Locus ID:	94071
Cytogenetics:	6 F3
Gene Summary:	Lectin-type cell surface receptor.[UniProtKB/Swiss-Prot Function]