

Product datasheet for MC206847

Ptgds (BC038083) Mouse Untagged Clone

Product data:

| | |
|---------------------------|--|
| Product Type: | Expression Plasmids |
| Product Name: | Ptgds (BC038083) Mouse Untagged Clone |
| Tag: | Tag Free |
| Symbol: | Ptgds |
| Synonyms: | PGD2, L-PGDS, 21kDa |
| Mammalian Cell Selection: | Neomycin |
| Vector: | PCMV6-Kan/Neo (PCMV6KN) |
| E. coli Selection: | Kanamycin (25 ug/mL) |
| Fully Sequenced ORF: | >MC206847 representing BC038083 Red=Cloning site Blue=ORF |

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCCGCATCGCC

ATGTGCAAGACAGTGGTAGCCCCCTCCACAGAAGGCGGCCTCAATCTCACCTCTACCTTCTCAGGAAAA
ACCAAGTGTGAGACCAAGATCATGGTACTGCAGCCTGCGGGGGCTCCTGGACACTACACCTACAGCAGCCC
CCACTCGGGCAGCATCCACTCCGTGTAGTGGTGGAGGCCAACTATGACGAGTACGCTCTGCTATTCAGC
AGAGGCACCAAGGGCCAGGCCAGGACTTCCGCATGGCCACCCTCTACAGCAGAACCAGACTCTGAAGG
ACGAGCTGAAGGAGAAATTCACCACCTTTAGCAAGGCCAGGGCCTCACAGAGGAGGACATTGTTTTCTT
GCCCCAACCGGATAAGTGCATTCAAGAGTAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Chromatograms: https://cdn.origene.com/chromatograms/ja1819_d01.zip

Restriction Sites: SgfI-MluI

ACCN: BC038083

Insert Size: 381 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

| | |
|-------------------------------|---|
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| Note: | Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required. |
| RefSeq: | <u>BC038083</u> , <u>AAH38083</u> |
| RefSeq Size: | 1354 bp |
| RefSeq ORF: | 381 bp |
| Locus ID: | 19215 |
| Cytogenetics: | 2 17.28 cM |
| Gene Summary: | Catalyzes the conversion of PGH2 to PGD2, a prostaglandin involved in smooth muscle contraction/relaxation and a potent inhibitor of platelet aggregation. Involved in a variety of CNS functions, such as sedation, NREM sleep and PGE2-induced allodynia, and may have an anti-apoptotic role in oligodendrocytes. Binds small non-substrate lipophilic molecules, including biliverdin, bilirubin, retinal, retinoic acid and thyroid hormone, and may act as a scavenger for harmful hydrophobic molecules and as a secretory retinoid and thyroid hormone transporter. Possibly involved in development and maintenance of the blood-brain, blood-retina, blood-aqueous humor and blood-testis barrier. It is likely to play important roles in both maturation and maintenance of the central nervous system and male reproductive system.[UniProtKB/Swiss-Prot Function] |