

Product datasheet for **MC206870**

Mrpl17 (BC010718) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Mrpl17 (BC010718) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Mrpl17
Synonyms:	MRP-L26; Rpml26
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC206870 representing BC010718. Blue=ORF Red=Cloning site Green=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAGGTTGTCGCTGGCTGCCGCTATCTCCACGGCCGCGTCTACCGCCGCTGGGCTGGGTCCGGAG
TCCCGCATCCATCTGCTGCGGAACCTTGCTCACGGGCCTAGTTCGCCACGAACGCATCGAGGCGACATGG
GCACGCGCGGACGAGATGAGGGGCTACGCGGAGAAGCTCATCGATTATGGAAAGCTGGGCGACACCAAC
GAACGAGCCATGCGTATGGCTGACTTCTGGCTCACAGTGAGTGCTCCCTACCTGCTAGTGAGGAAGGAA
TCCTTAGAGAACCGGTACCTGGCAGATTGAAGGCAGGGCATTCTTTTATGTACCAGACTTGCGAACTA
CCTTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-MluI



[View online »](#)

Plasmid Map:



ACCN: BC010718

Insert Size: 351 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC010718](#)

RefSeq Size: 1039 bp
RefSeq ORF: 350 bp
Locus ID: 27397
Cytogenetics: 7 E3
MW: 13.3 kDa