

Product datasheet for MC206878

Hebp2 (BC006898) Mouse Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	Hebp2 (BC006898) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Hebp2
Synonyms:	SOUL
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC206878 representing BC006898. Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCAGAGGAGCCGGAGCCAGATCTCGGGGTGGCCGAGGGCTCAGAGGATCAGGCCTTGGAGATGCCG
AGCTGGAAAGCTCCGGAGGACATAGACCCCGAGCAGGAAGTTATGAGATCCGGCACTACGGACCGCT
AAGTGGGTGAGCACTTGGTGGAGTCCCTGGACTGGGATTCAGCCATCCAGACTGGTTTACAAAAATTG
AATGGCTACATTCAGGCAAAAATGAGAAAGAGATGAAAATCAAGCTGACAGCCCGGTGACAAGCTAC
GTGGAGCCCGGCTCAAGTCCTTTTCAGTCTTTTGGATTCTCCAGTGGCCAAAAGAATCAAGAACAAC
TTTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: Sgfl-MluI



[View online »](#)

Plasmid Map:


ACCN: BC006898

Insert Size: 351 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in *E. coli* are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC006898](#)

RefSeq Size: 986 bp

RefSeq ORF: 350 bp

Locus ID: 56016

Cytogenetics: 10 A3

MW: 12.8 kDa

Gene Summary: Can promote mitochondrial permeability transition and facilitate necrotic cell death under different types of stress conditions (By similarity). May have low affinity for heme (PubMed:15518569).[UniProtKB/Swiss-Prot Function]