

## Product datasheet for **MC206906**

### 1300014I06Rik (BC037079) Mouse Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	1300014I06Rik (BC037079) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	1300014I06Rik
Synonyms:	1300014I06Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC206906 representing BC037079. Blue=ORF Red=Cloning site Green=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCCTCGGCGGTGTTTGAGGGCACGTCGCTCGTGAACATGTTTCGTGCGCGTTGCTGGGTGAACGGG
ATCCGCAGGTCATCGTCAGCAGGCGCGCGATGAGGAGGAGTTCTCGAAATCCGCACCGAGTGGTCG
GACCGCAGCGTCTACCTGCACCGCAGCTTCGCCGACCTGGGCGCCTGTGCAAGCGCCTCCGCGAT
GCCTTCCCGAGGACCGGTCTGAGCTGGCGCGCACGCCGCTGCGGCAAGGACTGCTCGCCATCAAGGT
GCACATGACATAGAGACCAGGCTCAATGAGGTCGAGAAGCTGCTGGAGACGGTCATCAGCATGCCGTGT
AAGTATCCAGGTCTGAGTTGTTCTCACCTTCTTTGAAAGATCGCCTTTGGATCAGGTGTTAAAAAT
GATAATGTCCATAAGATCCAACCCAGCTTCAAAGTCCTGTTAAATATCAGAGGCTAAAGAAGGCACT
GGAGCTCTCAGAACTGGAGTGGTAGACGGTTGTGAGCCACCAGGTGGGTGCTGGGAACTGAGTGTGGT
CTTCTGGAAGAACATCAAGTATTCTTAAGTCTGAGCGTATCTTCAGCCCCAGGAAAGCATTTCTGAGA
AATCATGAGGTCCAATGGATTTTGTGGCAAACACAGAGACCATAGTCATTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: Sgfl-MluI



[View online »](#)

## Plasmid Map:



ACCN: BC037079

Insert Size: 675 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC037079](#)

RefSeq Size: 2041 bp

RefSeq ORF: 674 bp

Locus ID: 66895  
Cytogenetics: 13  
MW: 25.7 kDa