

## Product datasheet for MC206922

### 0610037D15Rik (BC022733) Mouse Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	0610037D15Rik (BC022733) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	0610037D15Rik
Synonyms:	0610037D15Rik; 4933430J04Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC206922 representing BC022733. Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAACAGGAGCCCTAGCTGGTCTGATCAGCTTCATGCGGTTCTCAGTTCAACTGATGAAAAATGGCC
CGGATTAAGAGTTTGTATCCCTTTTCAGTCTTCTACAGGGGACCTGGCTGACACCTGGACTTCTCAT
CGCTTACCTCCTCAGCCAGAGGCCCAAGCTCAACAGCCTTGGGCTCTAAAGACTCCCCTTTGGAGAG
AGACTGAGCTGTACTGAAGCCCACAGCTCCTCTCTCTGGGATGAAGTTACTAACTCCGATCCCAGCTT
CAGGCTCAGGCTCAGGTACTGAGGCGCTGAAGCATTCTGTACAGAGCCTGTTGAAGGATCGGGAGCAA
CATATGTTCAAGATCAGTCTCTAGAAGCATCACTAAAAGTCTGCAGGAGGGCCCAAAGGAGAGCC
TTCTGGAGCAATGCCTGGAAGGCTGAAAAGGAACTGCAGGGCCTTCGGAGTCAGTTCAAGAGCTA
GTCCAAAACCAAATGAAAGCAAGACCAGGAAAGTTGGGTCTGCCAGTGGCTTTCACCAAGAGCTGCAG
AGTGAAATGACTGCTGTCACTCATGTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites:	Sgfl-MluI
ACCN:	BC022733
Insert Size:	582 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [BC022733](#)

**RefSeq Size:** 1251 bp

**RefSeq ORF:** 581 bp

**Locus ID:** 68394

**Cytogenetics:** 4 D1

**MW:** 21.6 kDa