

Product datasheet for MC206932

2310061I04Rik (BC100355) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: 2310061I04Rik (BC100355) Mouse Untagged Clone
Tag: Tag Free
Symbol: 2310061I04Rik
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >MC206932 representing BC100355.
 Blue=ORF Red=Cloning site Green=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGAGGAACACCTGGCTGTCATGCATGAGAGACTAAGACAAGAGCTTCCCACCCTTCTCCTTCGCTCC
CATGACTACACTATATATTCAATGGATGTGGAATTCATCAATGAGATCCTGAACATTCGTACCAAGGGC
CGGACATTTTATGTTATGTCCCTGACCCTTGGCCGTTCTGTTTGGAACTATTTGCACAGTTTCGG
TTGGAGATTTTACAGCTGACCCGCCACCCTGAGAACTGGACCCTGCAAGCCCGCTGGCGTCTAATAGGT
CTGCCATACATATGCTCTTCTGCGTTTCTACAGGCGTGATAAAGAGGATCTTTACCGGACCTTCGAT
GCCTTTCCACCTTCTACCTGAATTCAGTGGCCTCATCTGCGCCATCACCTGGACAAGCTGATGCCT
TCACACTACCGTCGACGCTGTGAAAAAGCTGCTGGTGGGAGCCCTGGTGGCCCTGGGGCTGTCAAGAA
CCAGAGCCAGCTTACACCTGTGCTCCAAGGCCTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:



ACCN: BC100355

Insert Size: 519 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC100355](#)

RefSeq Size: 893 bp
RefSeq ORF: 518 bp
Locus ID: 69662
Cytogenetics: 17 B1
MW: 20.5 kDa