

## Product datasheet for MC206949

### 1700052K11Rik (BC099431) Mouse Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	1700052K11Rik (BC099431) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	1700052K11Rik
Synonyms:	A1840954
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC206949 representing BC099431. Blue=ORF Red=Cloning site Green=Tag(s)


```
GTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGATCCGAGGAGAGAGCCCGGGAGCGGGCAAGGGGTGTCCTGACTGAGTTCGCGGTGCCACCCGATGT
CCGCCAGCCGGCCATCCCCGCAGCCATGGTCCACGTCCCAGACCCACGCCACGCCACGCCACGCCA
CAGGCACGAAAGCCGGCCGGCGCTCTCGCCTCCACGTACCTGACTGGCCTGGCGTCTCCAGGGACA
AAGGTGTCCTGGGATCACAGCCGAGCCCTACCAGCATGCCCCCTCCAGAGCGGAGGTCGCTACCCGC
CCCCGGGGAACCTCTCAAAGCACCTAGGGGGAGGTTCCATCTTAACCTCCGGTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

## Plasmid Map:



**ACCN:** BC099431

**Insert Size:** 333 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: BC099431  
RefSeq Size: 1329 bp  
RefSeq ORF: 332 bp  
Locus ID: 73431  
Cytogenetics: 11 E1  
MW: 11.8 kDa