

Product datasheet for **MC206977**

Fmr1nb (NM_001166619) Mouse Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	Fmr1nb (NM_001166619) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Fmr1nb
Synonyms:	3830422N12Rik; NY-SAR-35; NYSAR35
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC206977 representing NM_001166619. Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCCTTCGGATCGAAGGCCAAGTCAGAGAAGGAATAGATCTAAGAGCCGTGATTATCGTGGTGCACGG
TCAAAGGTAACAAGAGCTGATACGAGGAACAGAGACGATACTCTTGCCCTCAGTATGTATCAGGGCCCT
CCGAGTGCCGACCAGGGGAACAACATGGCGGATGCCCTCGGTTTGGCTTCTGGACTTCAGTAAGCCAG
TGTCTGCAATACTTGTGGCCAGGAGGCACTTGGCCCTGCTTCTACTTTTATTCTGGACGCTGGTGATC
CTGTTCCGTCCTGTGAACACTGCGAAATGCCCATTCTTGCTGAAGCTGCAGAACTGAACCCCTTTG
GGAAATATGTTGGACTTTTTCTTTCCAACAAGCCTACACAGGTGCTACGGGTGTTTGGCCTTGCTGCG
ATCAGCATTCTAGTCCTGGGATTCTGCCTATGTGCTGCTGCTCCATGTGCTGGAGGAGGAAGAGGATG
AACAGGATGTTGAAGTTTTGAAGAAACAGAAATCAAAGGGAAGAAGCCTAAAGGAAGGAAGGCGTCA
GAAGAGAGAGCTTTACTGTCCCATTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_001166619

Insert Size: 579 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001166619.1](#)

RefSeq Size: 820 bp

RefSeq ORF: 579 bp

Locus ID: 207854
UniProt ID: [Q80ZA7](#)
Cytogenetics: X A7.1
MW: 22 kDa