

## Product datasheet for MC206993

### Ccdc124 (BC050938) Mouse Untagged Clone

#### Product data:


Product Type:	Expression Plasmids
Product Name:	Ccdc124 (BC050938) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Ccdc124
Synonyms:	1810023B24Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC206993 representing BC050938. Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCCAAGAAGTTCAGGGCGAGAACAGCAAGTCGGCAGCTGCCGAGCACGGAGGGCTGAAGCCAAG
GCGGCCGCTGATGCCAAGAAACAGAAGGAGCTGGAAGATGCTTACTGGAAGGATGAGGACAAACACGTG
ATGCGGAAGGAGCAACGCAAGGAGGAGAAGGAGAAGCGACGCCTGGAGCAGCTGGAGCGCAAGAAGGAG
ACGCAGCGCCTGCTGGAGGAGGAAGACTCACGGCTCAAAGGCGGCAAGGCCCCACGTGTGGCGCCCGCC
AAGGTCACACGTGCACAGATCGAGGACTCGCTGCGGAGAGAGCAGCGGGCAGAGCCAGTGGAAAAGGCT
AAGAGCCACCTAGAGCTGCCACTGGAAGAGAACCTTAACCGGCGCCTGCAGGAAGAGGGCAGCGTGGAG
GCGCGCACTGTGGAGGACGCCATCGCAGTGCTCAGTGTGGCTGAGGAAGCTGACCGGCATCCTGAGCGC
CGGATGCGAGCGCCTTACC CGCTTTGAGGAGGTGCAACTGCCACGGTTGAAGCAGGAGAACCCCAAC
ATGCGGCTGTGCAACTGAAGCAGCTGCTGAAGAAGGAGTGGTTGCGCTCCCTGACAACCCGATGAAC
CAGCGGCCCTGCCCTTTAATGCGCCCAAGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: Sgfl-MluI



[View online »](#)

**Plasmid Map:**


**ACCN:** BC050938

**Insert Size:** 654 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [BC050938](#)

**RefSeq Size:** 912 bp

**RefSeq ORF:** 653 bp

**Locus ID:** 234388  
**Cytogenetics:** 8 B3.3  
**MW:** 25.3 kDa  
**Gene Summary:** Required for proper progression of late cytokinetic stages.[UniProtKB/Swiss-Prot Function]