

## Product datasheet for MC206994

### Tbpl1 (BC033926) Mouse Untagged Clone

#### Product data:


Product Type:	Expression Plasmids
Product Name:	Tbpl1 (BC033926) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Tbpl1
Synonyms:	4732475G08, TLF, TRP, TRF2, STD
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC206994 representing BC033926. Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGATGCAGACAGTGATGTTGCATTGGACATTTAATTACAAATGTAGTCTGTGTTTTAGAACAAGA
TGCCATTTGAACTTAAGGAAGATTGCTTTGGAGGGAGCAAATGTAATTTATAAGCGTGATGTTGGGAAA
GTATTAATGAAGCTTAGAAAACCTAGAATTACAGCTACAATTTGGTCCTCAGGAAAAATTATTTGCACT
GGAGCAACAAGTGAAGAAGAAGCTAAATTTGGTGCCAGACGTTTAGCCCGTAGTCTGCAGAACTAGGT
TTTCAGGTCATCTTCACAGATTTAAGGTGGTGAATGTTTGGCAGTGTGAACATGCCCTTTGAGATC
CGTTTGCCAGAATTTACAAAGAACAATAGACCTCATGCCAGCTATGAACCTGAACCTCATCCTGCCGTG
TGCTATCGGATAAAGTCTCTAAGAGCTACATTACAGATATTTTCAACAGGAAGCATCACAGTGACAGAT
CATGTAATGATCGAGTCCGTGGGCTAATTATTACTACATTGCATGGCCTCATCACCTGATGCCTGCC
CCAGAATCCTGCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-MluI



[View online »](#)

**Plasmid Map:**


**ACCN:** BC033926

**Insert Size:** 567 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [BC033926](#)

**RefSeq Size:** 1338 bp

**RefSeq ORF:** 566 bp

**Locus ID:** 237336

**Cytogenetics:** 10 A3

**MW:** 21.4 kDa

**Gene Summary:** Part of a specialized transcription system that mediates the transcription of most ribosomal proteins through the 5'-TCT-3' motif which is a core promoter element at these genes (By similarity). Seems to also mediate the transcription of NF1. Does not bind the TATA box. [UniProtKB/Swiss-Prot Function]