

Product datasheet for MC206995

Rpl12 (BC081469) Mouse Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	Rpl12 (BC081469) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Rpl12
Synonyms:	E430018F03
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC206995 representing BC081469. Blue=ORF Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCCGCCAAGTTCGACCCCAACGAGGTCAAAGTCGTGTACCTGAGGTGCACCGGAGGCGAGGTCGGC
GCCACATCCGCCTTGCCCCGAAGATCGGTCTCTGGGTCTGTCTCCGAAGAAAGTTGGCGATGACATT
GCCAAGGCTACCGGTGACTGGAAGGGTCTCAGAATTACAGTGAACTGACCATCCAGAACAGACAGGCC
CAGATCGAGGTGGTGCCTTCTGCCTCAGCCCTGATCATCAAAGCCCTCAAGGAGCCACCGAGAGACAGG
AAGAAGCAGAAGAACATTAACACAGTGGGAACATTACTTTTGATGAGATTGTCAACATTGCCCGGCAG
ATGAGACACCGGTCTTTGGCTAGAGAACTTTCTGGAATCAAGGAGATCCTGGGTACTGCACAGTCT
GTGGGCTGCAATGTGGATGGCCGCCACCCTCATGACATCATAGATGACATCAACAGTGGTGCAGTGGAG
TGCCAGCTAGTTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC

Restriction Sites: Sgfl-Mlul



[View online »](#)

Plasmid Map:


ACCN: BC081469

Insert Size: 498 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC081469](#)

RefSeq Size: 640 bp

RefSeq ORF: 497 bp
Locus ID: 269261
Cytogenetics: 2 B
MW: 17.8 kDa
Gene Summary: Binds directly to 26S ribosomal RNA.[UniProtKB/Swiss-Prot Function]