

Product datasheet for **MC207062**

Npm3 (BC048491) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Npm3 (BC048491) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Npm3
Synonyms:	Nub1
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
Fully Sequenced ORF:	>BC048491 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCGGCCGGCGCGCGGCTGCTTTAGCGTTCTTGAACCAAGAGAGCCGAGCCGAGCTGGGGCGTCC
GGGGCCTGAGGGTCCCGGCCCGGTCAATGGACAGTTCTTCTCGTTGTGAGCTCTCCGGCCACAC
CCGCTCCTTCAAGTTGAGGAGGAGGATGACACGAGCACGTGCTGGCGTTGAACATGCTGTGC
TTACCGAGGGGGCCACGACGAGTGAATGTGGTGGAAGTTGTGGCCCGGGATCACGACAATCAGGAGA
TTGCAGTTCTGTGGCAACCTCAGTTGTCTTCCAGCCATGCTCAGTGTGGATGATTTCCAACCTCCA
ACCACCTGTCACCTTCCGCCTGAAGTCCGGTCTGGCCAGTTCGGATCACTGGCCGGCATCAGATCGTT
TGTATAAACAACGACCTTCTGAAGAAGAAAGTATGATGAGAGCGAAGAAGATGAAATTAAGCTGTGTG
GCATCCTTCTGCCAAGAAGCACAGGGGCAGGCCTTAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCTGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	BC048491
Insert Size:	528 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>BC048491</u> , <u>AAH48491</u>
RefSeq Size:	793 bp
RefSeq ORF:	527 bp
Locus ID:	18150
Cytogenetics:	19 38.75 cM
Gene Summary:	Plays a role in the regulation of diverse cellular processes such as ribosome biogenesis, chromatin remodeling or protein chaperoning. Modulates the histone chaperone function and the RNA-binding activity of nucleolar phosphoprotein B23/NPM. Efficiently mediates chromatin remodeling when included in a pentamer containing NPM3 and NPM. [UniProtKB/Swiss-Prot Function]