

Product datasheet for **MC207088**

1810029B16Rik (BC016246) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: 1810029B16Rik (BC016246) Mouse Untagged Clone
Tag: Tag Free
Symbol: 1810029B16Rik
Synonyms: 1810029B16Rik
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >BC016246
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAAAAACAAGATATTCAAAGAAAGATGCCTGTGAATTAGTTGAAAGGTAAGTGGATCGGTTTCAGCA
GTGAGCTGGAGCAGATTGAGTTACAGAACAGCATCAAAGACAGGCAGGGAAGCGGCACCATTCCCGGGA
GGCCGTCAAGCAGACCCTGGAGCGGGAGCGGCAGCAGTACGAAGGCTATGGCTTTGAGATCCCAGAC
ATTTTAGACTCTAATACTCTGCAGACTTCCCGGAATGGGATTTTGATCTGAAGAAATTGCCAAACATCA
AAATGAGGAAACTTTGTGCTGATGATGCAGTTCCTAAGAAGCGCAACAGAAAAATATTTTAAATATAGA
AAAGGATTTAGGAGAATTAGAACTATCAGGGCCAACTGGAGCCACTACAGATGGGAACTGGAACCGCA
AGTGAATCAAGTGACACGGATGAAGAAATGACTCCAGTGCCTGTCAGCCCTCATTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCTGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-MluI
ACCN: BC016246
Insert Size: 477 bp
OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC016246](#), [AAH16246](#)

RefSeq Size: 3162 bp

RefSeq ORF: 476 bp

Locus ID: 66282

Cytogenetics: 8 B3.1