

Product datasheet for MC207162

Tpbpb (BC099441) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Tpbpb (BC099441) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Tpbpb
Synonyms:	1600012N09Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC099441

Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGCATCGCC

ATGACTCCTACAGTCTTCTTGTCATCCTGTGCCTGGGTGTGGCCTCAGCTGCTATAATCCCTGAACCAA
 CTTTGGATACTGAAGTGAAGAGCAGAAGGATAAAGAAGGACTCAAAAAAGCAACGTGGAATGAGTTTGT
 CATGAACTGAACAATAGCAAAATGACCAGGAGAAGGATGGCATCGACATAGAATTGAGAGTCTTCTTT
 GGTCAGCTGACTGATGAAGAACTAATGAAAATAATGACTAACATTTTCCACCAATGTTGGAAGAGGAGA
 AAACACAGCCAGTTGTTGATGACCCTGAATTTGAGGATTACACAGAGAGTGGCGATGGGTTTTTTGTACC
 AAATCAGGTATGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCTGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	BC099441
Insert Size:	774 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>BC099441</u>
RefSeq Size:	1017 bp
RefSeq ORF:	362 bp
Locus ID:	116913
Cytogenetics:	13 B2