

## Product datasheet for **MC207188**

### Serpina4 (BC031891) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Serpina4 (BC031891) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Serpina4  
**Synonyms:** KALps, Klst  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Fully Sequenced ORF:** >BC031891

Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGCTGTCCTTAGAGGCTGTCGGGGACACCCAGACTCAGATCCTCGAAGGCCTGGGCTTCAACCTCACAA  
AGCTGTCAGTGCCTGAAATCTATGAGGGTTTCAGGTCCATTCAAGACACGTTAGCTCGCCATTCACTGA  
GCTGAACATCTCTGTAGGCAATGCACTGATCCTGAGCCAGGACCTGCGACTCCTCCATGAGTTCGCCAGT  
CCCACCGAATCCTCCTTTAACAGCAAAGTCTTACAAGCTAACTTCGGAGACGAGAAGGCTGCGATCCAGC  
TCATCAACGACTATGTCAAGGAGAAAAGTCAGGGGAAGATCAAAGATTTGATCAGTAAACTGAGCCAAA  
TGTAAGGATGGTCTGGTGAATTATATCTTCTTCCAAGGCCTTTGGGAGAAACCCTTACATACTTCTGGA  
ATCTCTACCCAAGACTTCTACGTAGATGAGAACACAGTTGTGAAGGTACCAATGATGAGGCAAGTTCAGA  
AACATGGGTATCTTGAGGACAGCAATGTACCCTGCAAAGTCTGCGGATGGATTACAGGGGAGATGCAAT  
GGCGTTCTTATCCTTCTAACCCTGGGCAAGATGAATGAGGTGGAGAAGGTGCTGTCCCCAGGCATGCTA  
ACAAAGTGGCATCGCTTGCTTCAAGATAGGTAA

**ACGCGT**ACGCGGGCCGCTCGAGCAGAAACTCATCTCTGAAGAGGATCTGGCAGCAAATGATATCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** Sgfl-MluI  
**ACCN:** BC031891  
**Insert Size:** 636 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">BC031891</a></u>
<b>RefSeq Size:</b>	3933 bp
<b>RefSeq ORF:</b>	662 bp
<b>Locus ID:</b>	321018
<b>Cytogenetics:</b>	12 E