

Product datasheet for MC207234

Fxyd2 (NM_007503) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Fxyd2 (NM_007503) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Fxyd2
Synonyms:	Atp1g1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC207234 representing NM_007503 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCTGGGAAATATCAGATCTGTCAGCGAACAGTGGTGGCAGTGCCAAGGGGACAGAGAATCCCTTCG
AGTACGACTATGAAACCGTCCGCAAAGGAGGCCTGATCTTCGCGGGCCTGGCCTTCGTCGTTGGCCTCCT
CATCATTCTCAGCAAAGGTTCCGCTGTGGGGCGGTAAGAAACATAGGCAGGTCAATGAAGATGAAGT
TGA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_007503
Insert Size:	213 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_007503.3](#), [NP_031529.2](#)

RefSeq Size: 569 bp

RefSeq ORF: 213 bp

Locus ID: 11936

UniProt ID: [Q04646](#)

Cytogenetics: 9 A5.2

Gene Summary: This gene encodes a member of a family of small membrane proteins that share a 35-amino acid signature sequence domain, beginning with the sequence PFXYD and containing 7 invariant and 6 highly conserved amino acids. The approved human gene nomenclature for the family is FXYD-domain containing ion transport regulator. Mouse FXYD5 has been termed RIC (Related to Ion Channel). FXYD2, also known as the gamma subunit of the Na,K-ATPase, regulates the properties of that enzyme. FXYD1 (phospholemman), FXYD2 (gamma), FXYD3 (MAT-8), FXYD4 (CHIF), and FXYD5 (RIC) have been shown to induce channel activity in experimental expression systems. Transmembrane topology has been established for two family members (FXYD1 and FXYD2), with the N-terminus extracellular and the C-terminus on the cytoplasmic side of the membrane. The Type III integral membrane protein encoded by this gene is the gamma subunit of the Na,K-ATPase present on the plasma membrane. Although the Na,K-ATPase does not depend on the gamma subunit to be functional, it is thought that the gamma subunit modulates the enzyme's activity by inducing ion channel activity. Multiple transcript variants have been described for this gene that are expressed in tissue-specific and developmental stage-specific patterns and encode proteins that differ at the N-terminus. [provided by RefSeq, Sep 2009]

Transcript Variant: This variant (a) represents the most frequently occurring transcript and it encodes the longest protein (isoform a).