

Product datasheet for MC207266

Aimp1 (NM_007926) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Aimp1 (NM_007926) Mouse Untagged Clone
Tag: Tag Free
Symbol: Aimp1
Synonyms: 9830137A06Rik; AIMP1/p43; Emap2; EMAP11; p43; Scye1
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC207266 representing NM_007926
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGC**C

ATGTTTCTCTGCCGATTCTGGGAAAGATGGCAACGAATGATGCTGTTCTGAAGAGGCTGGAGCAGAAGG
 GTGCAGAGGCGGATCAGATCATGAATATCTCAAGCAGCAGGTTGCTCTTCTTAAGGAGAAAGCAATTTT
 GCAGGCAACAATGAGAGAAGAAAAGAACTTCGAGTTGAAAATGCTAAACTGAAAAAGAAATAGAAGAG
 CTAAGGCAAGAGCTGATTCTGGCAGAAATTCATAACGGAGTGGAGCAAGTGCCTGTTTCGATTGAGTACTC
 CACTGCAGACGAACTGTACTGCTTCTGAAAAGTGTGGTGCAGTCTCCATCAGTAGCAACCACCGCCTCTCC
 TGCTACAAAAGAGCAGATCAAAGCGGGAGAAGAAAAGAGGTGAAAGAGAAGACTGAAAAGAAAGGAGAG
 AAAAGGAGAAGCAGCAGTCCGACGAGCAAGTACTGACTCCAAGCCTATCGACGCATCGCGTCTGGATC
 TTCGAATTGGTTGTATTGTTACTGCCAAGAAGCACCTGATGCAGATCACTGTATGTGGAGGAAGTAGA
 TGTGGGAGAAGCAGCCCCGCGCACGGTCGTCAGCGGGCTGGTGAATCATGTTCTCTAGAACAGATGCAA
 AATCGTATGGTGGTTTTACTCTGTAATCTGAAGCCTGCAAAGATGCGGGGAGTTCTGTCTCAAGCCATGG
 TGATGTGTGCCAGTCCACAGAAAAGTGGAGATTCTGGCCCTCCCAACGGGTCGGTCTCTGGGACAG
 AATTACTTTTGTGCTTTTCTGGAGAGCCTGACAAGGAGCTAAACCCTAAGAAGAAGATCTGGGAGCAG
 ATCCAGCCTGACCTGCACACCAATGCTGAGTGTGTGGCCACATACAAAGGAGCTCCCTTTGAGGTGAAGG
 GGAAGGGAGTTTGCAGAGCCCAAACCATGGCCAATAGTGAATTAATA**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_007926
Insert Size: 960 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_007926.2 , NP_031952.2
RefSeq Size:	1122 bp
RefSeq ORF:	960 bp
Locus ID:	13722
Cytogenetics:	3 G3