

## Product datasheet for MC207320

### Lyl1 (NM\_008535) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Lyl1 (NM\_008535) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Lyl1  
**Synonyms:** bHLHa18; Lyl-1  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC207320 representing NM\_008535  
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGTGCCCCCCCCAGGCCAGAGCAGAGGTGGGTTCCGCCATGACTGAGAAAAGTGGATGGTATGTGCCT  
 CCAGCCCAGCACCTGCCCTCCCTCTAAGCCGGCCTCACCTGGGCCCTTTCTACAGAGGAGGTGGACCA  
 CCGAAACAGTGCCTCCCTGGTTGCCTCCCGGGTGGCAGTATAAACCTGGGACACACCAGGCCATA  
 GGGGCAGCCATGCCACCACAGAGCTCAGTGCCTTTCCGGCCCTCCTTGCTGCAGCTAACTGCCTTGGAA  
 GAGCTCCACCTACCCTGGCTGTACATTACCACCCTCACCCCTTCTCAACAGTGTCTACATTGGGCCAGC  
 AGGACCCTTCAGCATCTCCCTAACAGCCGGCTGAAGCGCAGACCAAGCCATAGTGAGCTGGACTTGGCT  
 GACGGGCACCAACCCAGAGGTGGCTCGGCGAGTGTTCACCAACAGCCGTGAGCGCTGGCGGCAGCAGC  
 ACGTTAACGGCGCCTTCGAGAGCTCAGGAAGTGTGCCACCCACCCGCCGACCGGAAGTGAAGCA  
 GAACGAGGTGCTGCGCTGGCCATGAAGTATATAGGCTTCTGGTGGCTGCTGCGAGACCAACCGGT  
 GTGCTGACCTCCGGCCCCAGCGCTCCCGGTCCCGCAAGCCACCTGCGCGCAGGGCGTGGAGGGCAGCG  
 CACGCTTCGGGGCCGGGCACAGAGTAGAGGTGCACGCTCACAGCCGTGCTTCTGGGGACTGTGATGG  
 CGACCCCAATGGGTCAGTGAGACCCATCAAGTTGGAGCAGACGTCCCTGAGTCTGAGGTGCGGTGA

ACGCGTACGCGGCCGCTCGAGCAGAAAAGTCTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_008535  
**Insert Size:** 837 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_008535.2</a> , <a href="#">NP_032561.2</a>
<b>RefSeq Size:</b>	1818 bp
<b>RefSeq ORF:</b>	837 bp
<b>Locus ID:</b>	17095
<b>UniProt ID:</b>	<a href="#">P27792</a>
<b>Cytogenetics:</b>	8 41.02 cM