

## Product datasheet for **MC207361**

### **Psmb10 (NM\_013640) Mouse Untagged Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Psmb10 (NM\_013640) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Psmb10  
**Synonyms:** Mecl-1; Mecl1  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC207361 representing NM\_013640  
**Red=Cloning site Blue=ORF Orange=Stop codon**

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGCTGAAGCAGGCAGTGAACCCACAGGAGGCTTCTCTTTCGAGAACTGCCAGAGGAATGCGTCCTTGG  
AACACGTCCTTCCGGGACTTCGGGTTCCCTCATGCACGCAAGACCGGGACTACCATCGCGGGGCTTGTGTT  
CCGAGATGGAGTCATCCTGGGAGCGGACACGCGGCCACTAACGATTCGGTTGTGGCGGACAAAAGCTGC  
GAGAAGATCCACTTCATCGCCCTAAAACTACTGCTGTGGGGCTGGAGTAGCTGCGGACACTGAGATGA  
CTACGCGGATGGCAGCTTCCAAGATGGAACACATGCGCTGTCCACCGGCCGTGAGCCTCGGGTGGCCAC  
GGTCACCCGATCCTGCGCCAGACGCTTTTCCGGTACCAAGGCCACGTGGGGGCATCACTGGTCGTGGGC  
GGGGTTGATTTGAACGGACCTCAGCTCTACGAGGTGCACCCACATGGTTCCTACAGCCGTCTGCCCTTA  
CTGCCCTTGGCTCTGGTCAGGGTGCAGCCGTGGCACTGTTGGAAGACCGTTCCAGCCAAACATGACGCT  
GGAGGCTGCGCAAGAGCTGTTGGTGAAGCCATCACAGCAGGGATCCTGAGTGACCTGGGCTCTGGGGGC  
AATGTGGATGCCTGTGTGATCACTGCAGGGGTGCCAAGCTGCAGAGAGCATTGAGCACCCACCAGC  
CTGTGCAGAGAGCTGGCCGTTACCGATTTGCTCCTGGAACCACACTGTCTGACCCGGGAAGTGAGACC  
CCTGACCCTGGAATCTTGAGGAACTGTGCAGGCCATGGAGGTGGAATGA

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_013640  
**Insert Size:** 822 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_013640.3</a> , <a href="#">NP_038668.2</a>
<b>RefSeq Size:</b>	1270 bp
<b>RefSeq ORF:</b>	822 bp
<b>Locus ID:</b>	19171
<b>UniProt ID:</b>	<a href="#">O35955</a>
<b>Cytogenetics:</b>	8 53.06 cM
<b>Gene Summary:</b>	The proteasome is a multicatalytic proteinase complex which is characterized by its ability to cleave peptides with Arg, Phe, Tyr, Leu, and Glu adjacent to the leaving group at neutral or slightly basic pH. The proteasome has an ATP-dependent proteolytic activity. This subunit is involved in antigen processing to generate class I binding peptides. Plays a role in determining the T-cell repertoire for an antiviral T-cell response.[UniProtKB/Swiss-Prot Function]