

Product datasheet for **MC207374**

Rad51d (NM_011235) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Rad51d (NM_011235) Mouse Untagged Clone
Tag: Tag Free
Symbol: Rad51d
Synonyms: R51H; R51H3; Rad5; Rad51I3; TRAD; Trad-d5
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC207374 representing NM_011235
Red=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGGCATGCTCAGGGCAGGGCTGTGCCCGGGCTCACCGAGGAGACCGTCCAGCTTCTCAGAGGCCGAA
AGATAAAACAGTGGCAGACCTGGCAGCTGCTGACTTGGAGGAAGTAGCCCAAGAAGTGTGGCTTGCCTA
CAAGGCCCTCGTTGCCCTGAGGAGGGTGTGCTGGCGCAGTTCTCGGCTTCCCCTTAAATGGCGCAGAT
CTCTATGAGGAACTGAAGACTTCCACGGCCATCCTGTCCACCGGCATCGGAAGCCTGGACAACTACTTG
ATGCTGGCCTCTATACTGGGAGGTGACTGAAATTGTGGGTGGCCAGGTAGCGGCAAAACCCAGGTGTG
TCTCTGTGTGGTGCAAAATGTGGCCATAGCCTGCAGCAGAATGTACTGTATGTGGATTCCAATGGAGGA
ATGACGGCGTCCCGCCTCCTCCAGCTACTACAGGCTAGAACCAGATGAGGAGAAAACAGGCAAGTGCTC
TCCAGAGGATACAGGTGGTGCCTTCATTTGACATCTTCCGGATGCTAGATATGCTACAGGACCTTCGCGG
CACCATAGCCCAGCAGGAAGCAACTTCTCAGGCGCCGTGAAGGTTGTGATTGTGGACTCGGTCACTGCA
GTGGTCCGCCCCACTTCTGGGAGGTCAGCAGAGGGAAGGCCTGGCCTTGATGATGCAGCTGGCCCCGAGAGC
TCAAGATCCTGGCCGGGACCTGGGTGTGGCAGTGGTGGTGACCAACCACTTGACTCGAGATTGGGATGG
TAGAAGATTCAAACCTGCCCTTGGAGCCTCCTGGAGCTTTGTGCCAGTACCCGGATTCTCCTGGATGTC
ACTGAGGGGGCTGGGACTCGGTAGCAGCCAACGCACAGTATGTCTGACCAAGTCTCCCCGCCAGCCAA
CGGGTCTGCAGGAGATGATAGACATTGGGACATTGGGACTGAGGAGCAGAGCCAGAATTACCTGGCAA
GCAGACGTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_011235



[View online >](#)

Insert Size:	990 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_011235.4 , NP_035365.1
RefSeq Size:	7072 bp
RefSeq ORF:	990 bp
Locus ID:	19364
UniProt ID:	O55230
Cytogenetics:	11 50.3 cM
Gene Summary:	<p>This gene belongs to the Rad51 gene family whose products play a major role in homologous recombination and DNA repair. The encoded protein interacts with other proteins of this family, including Rad51b, Rad51c and Xrcc2, and plays an essential role in both DNA repair and telomere maintenance. In humans, germline mutations in this gene may be associated with predisposition to ovarian cancer. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, May 2013]</p> <p>Transcript Variant: This variant (1, also known as Rad51d-FL) encodes the longest isoform (1).</p> <p>Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>