

## Product datasheet for **MC207385**

### **Rps7 (NM\_011300) Mouse Untagged Clone**

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Rps7 (NM_011300) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Rps7
Synonyms:	Mtu; Rps7A; S7
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
Fully Sequenced ORF:	>MC207385 representing NM_011300 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGTTTCAGCTCGAGCGCCAAGATCGTGAAGCCCAATGGCGAGAAGCCGGACGAGTTCGAGTCTGGCATCT  
CTCAGGCGCTGCTCGAGCTGGAGATGAACTCCGATCTCAAGGCGCAGCTGCGGAACTCAACATCACCGC  
AGCCAAGGAAATCGAAGTTGGTGGTGGTCGAAAGCTATCATAATTTTTGTACCAGTTCCTCAGCTTAA  
TCTTTCCAGAAAATCCAAGTCCGGCTAGTTCGTGAATTGGAGAAAAAGTTCAGTGGAAGCACGTTGGTCT  
TCATTGCTCAGAGGAGGATTCTGCCAAGCCAACCCGAAAAAGCCGTACGAAAAATAAGCAGAAGCGCCC  
CAGAAGCCGTACCTGACAGCAGTGCACGACGCCATCCTTGAGGACTTAGTCTTCCCAAGTAAAATTGTG  
GGTAAGAGGATCCGTGTGAACTGGATGGCAGCCGGCTCATAAAAGTTCACCTTAGACAAAGCACAGCAGA  
ACAACGTGGAACACAAGTTGAACTTTTTCTGGTGTGTACAAGAAGCTCACAGGCAAGGATGTTAATTT  
TGAATCCAGAGTTTCAGTT**TAA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_011300
Insert Size:	585 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_011300.3</a></u> , <u><a href="#">NP_035430.1</a></u>
<b>RefSeq Size:</b>	961 bp
<b>RefSeq ORF:</b>	585 bp
<b>Locus ID:</b>	20115
<b>UniProt ID:</b>	<u><a href="#">P62082</a></u>
<b>Cytogenetics:</b>	12 A2
<b>Gene Summary:</b>	Required for rRNA maturation.[UniProtKB/Swiss-Prot Function]