

Product datasheet for MC207392

Sdf4 (NM_011341) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Sdf4 (NM_011341) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Sdf4
Synonyms:	Cab4; Cab45
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
Fully Sequenced ORF:	>MC207392 representing NM_011341 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGTCTGGCTGGTGGCAATGACGCCAGGCAGAGTTCCTCTGTGGTCTAGCTGCTCATGGCCTCTGGT
TCTTGGCCCTTGCTCTTCTGATGGATGCAACCGCTAGACCTGCCAACCACTCATCTACTCGGAAAGAGC
GGCTAACAGGGAGGAAAATGAGATCATGCCCCAGACCACCTGAATGGGGTGAAGCTGGAGATGGATGGA
CACCTCAACAAGGACTTCCATCAGGAGTTTTCTGGGAAAGGACATGGATGGGTTTGATGAAGACTCAG
AGCCACGGAGAAGCCGGAGGAAGCTGATGGTCATCTTTCCAAGGTAGACGTGAACACTGACCGGAGGAT
CAGCGCTAAGGAGATGCAGCACTGGATTATGGAGAAAACAGCAGAGCACTCCAGGAGGCCGTCAAGGAA
AACAACTGCACTTCAGGGCTGTGGACCCTGACGGTGACGGCCATGTTTCTGGGATGAATATAAAGTGA
AGTTTTTGGCAAGCAAAGGCCACAATGAAAGGGAGATTGCTGAAGCCATCAAGAACCATGAGGAGCTCAA
AGTGGATGAGGAGACACAGGAAGTCCTTGGGAACCTCAGAGCCGATGGTATCAGGCAGACAATCCTCCT
GCAGACCTGCTGCTGACTGAGGACGAGTTCCTTTCATTCCTTACCCTGAGCACAGCCGGGCATGCTCA
AGTTCATGGTCAAGGAGATCTTTCGGGACTTGGATCAGGATGGTGATAAGCAGTTATCTCTGCCCGAGTT
CATCTCTCTGCCTGTGGTACTGTTGAGAACCAACAAGGCCAAGACATTGATGACAACTGGGTGAAGAC
AGAAAGAAGGAGTTTGAAGAACTGATTGACTCTAACCATGATGGGATTGTGACCATGGAGGAGCTAGAGA
ACTACATGGACCCATGAATGAATAACAATGCCCTCAACGAAGCCAAACAGATGATTGCCATTGCTGATGA
GAACCAGAACCACCCTGGAGCCCGAGGAGATCCTCAAGTACAGTGAGTTCTTACCAGCAGCAAGCTC
ATGGACTATGCCCGCAACGTGCATGAAGAGTTC**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_011341



[View online >](#)

Insert Size:	1086 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_011341.5, NP_035471.1</u>
RefSeq Size:	5255 bp
RefSeq ORF:	1086 bp
Locus ID:	20318
UniProt ID:	<u>Q61112</u>
Cytogenetics:	4 E2
Gene Summary:	<p>This gene encodes a member of the CREC family. The encoded protein contains multiple calcium-binding EF-hand motifs. This protein localizes to the Golgi lumen and may be involved in regulating calcium dependent cellular activities. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Oct 2014]</p> <p>Transcript Variant: This variant (1) represents the longest transcript and encodes the longer isoform (1). Variants 1, 2 and 3 encode the same protein.</p>