

Product datasheet for MC207402

Snrpd1 (NM_009226) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Snrpd1 (NM_009226) Mouse Untagged Clone
Tag: Tag Free
Symbol: Snrpd1
Synonyms: AA407109; AL023031; SMD1
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC207402 representing NM_009226
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGAAGCTCGTGAGATTTTTGATGAAATTGAGTCATGAACTGTAACCATTGAATTGAAGAATGGGACAC
 AAGTCCATGGAACAATCACAGGTGTAGATGTCAGCATGAACACACACCTTAAAGCTGTGAAAATGACCCT
 GAAGAACAGAGAACCCTGTACAGTTGAAACATTGAGTATTCGAGGCAATAACATTCGGTATTTTATTCTA
 CCAGACAGCTTACCTCTAGATACACTACTTGTGGATGTTGAACCTAAGGTGAAGTCTAAGAAAAGAGAAG
 CTGTTGCAGGAAGAGGCCGAGGCCGGGTAGAGGAAGAGGACGTGGTCTGGCAGAGGAAGAGGGGGTCC
 TAGGCGA**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-MluI
ACCN: NM_009226
Insert Size: 360 bp
OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online >](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_009226.4](#), [NP_033252.1](#)

RefSeq Size: 837 bp

RefSeq ORF: 360 bp

Locus ID: 20641

UniProt ID: [P62315](#)

Cytogenetics: 18 A1

Gene Summary: Plays role in pre-mRNA splicing as core component of the SMN-Sm complex that mediates spliceosomal snRNP assembly and as component of the spliceosomal U1, U2, U4 and U5 small nuclear ribonucleoproteins (snRNPs), the building blocks of the spliceosome. Component of both the pre-catalytic spliceosome B complex and activated spliceosome C complexes. Is also a component of the minor U12 spliceosome. May act as a charged protein scaffold to promote snRNP assembly or strengthen snRNP-snRNP interactions through non-specific electrostatic contacts with RNA.[UniProtKB/Swiss-Prot Function]