

Product datasheet for MC207405

Serpinb9c (NM_011453) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Tag Free
Symbol:	Serpinb9c
Synonyms:	3830421J05Rik; Nk9; Spill
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC207405 representing NM_011453 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCTGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGTTATGTAATAACAACCCCTTCCAAAAATGTATGCTATTCTCCCATCAACATCTCCTCTGCTCTAGCTA
TGTTCCCTTTGGGAGTGAAAGGAAACACCGAAATCCAGATATCTGAGGCAATTGGCTTAAACACAGCAAT
AGACATTATCAGAGCTTCCTGTGGATTCTCAATATCCTGAAGAAGCCAGCAAGTACACCTTCAGA
ATGGCCAAATAGGCTCTTTGCAGAGAACCTTGTGAATTCCTCCCAACATTTAAGGAGCCCTGTCTTCAGT
TCTATCACTGGGAGATGGAGCACCTCCCCTTACCAAGGCTCCAGAGGAGGCCAGGAATCACATAAACAC
ATGGGTCTGCAAAAATACAAAAGGGAATTCAGAAATTGTTGTCAAGTGGCTCGGTTGATTCAGAGACC
AGGCTTGTAATTGTCAATGCCTTATATTTCAAAGGAAGATGGCATCACCAGTTTGACATAAAGTCCACCA
GGAAAATGCCCTTTAAATAAACAAGGATGAGGAAAGACCAGTGCAGATGATGTTTCAGGAAGACATGTT
TAAGCTCGCCTATGTGAATGAGGTGCAGGTTCAAGTGTGTTGCTGCCCTATAAGGGCAAGGAGCTGAGC
TTGGTGGTCTGCTCCCGGATGATGGTGTGGAACCTCAGCAAGGTGGAAGCAATCTCACTTTTGAGAAGT
TATCAGCTTGACCAAACAGACTATTTGAAGACCACTAAAGTTTGGTTTTCCTTCCAAAATTTAACT
GGAAGATTATTATGACATGGAGTCCATCTTTCAGGACTTGGGAGTGGGGATATTTTCAAGGGGGCAAG
GCTGACTTATCAGAAATGTCTCCAGAGAGAGGCTGTGTGTGTCCAAGTTTCATTCAAAAGTGTGTTGTGG
AAGTTAATGAAGAAGGCACAGAAGCTACAGCAGCCACAGCTGATGATACTGTATGTTTCAGCTGAAACTCA
TGATGGCCAGACTTTCTGTGCTGACCACCTTTCTTTTCTTATCAGGCACAACAAACCAACAGCATC
CTGTTCTGTGGCAGTTCTCATTTCATAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



ACCN:	NM_011453
Insert Size:	1080 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_011453.6, NP_035583.1</u>
RefSeq Size:	1903 bp
RefSeq ORF:	1080 bp
Locus ID:	20707
Cytogenetics:	13 13.82 cM