

## Product datasheet for **MC207458**

### Spry1 (NM\_011896) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Spry1 (NM\_011896) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Spry1  
**Synonyms:** sprouty1; spry-1  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC207458 representing NM\_011896  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGGATCCCCAAGTCAGCATGGCAGCCACACTTCGCTAGTGGTGATTACGCCACCGGCTGTGGAAGGCC  
GGCAGAGGTTAGACTATGACAGGGACACTCAGCCTGCTACGATTCTGTCCTAGACCAGATCAAAGCCAT  
CAGAGGCAGCAATGAATACACAGAGGGACCTTCGGTAGCGAGAAGACCTGCTCCTCGCACTGCACCAAGA  
CCCGAAAAGCAGGAAAGGACTCATGAAATCATACCAGCCAATGTGAATAGCAGCTACGAGCACCGACCTG  
CCAGCCACCCGGGCAACGCCAGGGGCTCGGTGCTGAGCAGGTCCACCAGCACCGGAAGCGCGGCCAGCTC  
AGGGAGCAGCAGCAGTGTGTCTTCTGAGCAGGGCCTATTAGGACGGTCTCCGCCACCGGCCATCCCCA  
GGTCATAGGTCAGATCGGGTCAATCCGGACCCAGCCCAAGCAGCTGCTGGTGAAGACTTGAAGGCCTCCT  
TGAAAGAGGACCCACCCAGCACAAAGTTCATCTGCGAACAGTGTGGCAAGTGCAAATGTGGAGAGTGAC  
GGCCCCCGGGCCCTGCCCTCCTGTCTGGCCTGTGATCGGCAGTGCCTCTGCTCCGCGGAGAGCATGGTG  
GAATACGGGACCTGCATGTGCTGGTCAAGGGCATTCTTACCCTGCTCCAATGATGATGATGGAGGTT  
CTTACTCGGATAACCCATGCTCCTGTTCCAGTCCCAGTCCAGTCCAGATACCTGTGCATGGGAGCCCT  
GTCTTTGTGCCTACCCTGCTTGCTCTGCTACCCTCCCGCCAAGGGCTGCCTGAAGCTGTGCAGGGGCTGT  
TACGACTGGACCCACCGCCTGGCTGCAGGTGTAGAACTCCAACACTGTCTATTGTAAGCTGGAGAGCT  
GCCCTCAAGGGCTCAGGGCAAAGTGTCA**TGA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_011896  
**Insert Size:** 942 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_011896.3</a> , <a href="#">NP_036026.1</a>
<b>RefSeq Size:</b>	2493 bp
<b>RefSeq ORF:</b>	942 bp
<b>Locus ID:</b>	24063
<b>UniProt ID:</b>	<a href="#">Q9QXV9</a>
<b>Cytogenetics:</b>	3 B
<b>Gene Summary:</b>	May function as an antagonist of fibroblast growth factor (FGF) pathways and may negatively modulate respiratory organogenesis.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) represents the longest transcript. Variants 1, 2, 3, and 4, encode the same protein. Sequence Note:.