

Product datasheet for **MC207470**

Srp9 (NM_012058) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Srp9 (NM_012058) Mouse Untagged Clone
Tag: Tag Free
Symbol: Srp9
Synonyms: 9kDa
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC207470 representing NM_012058
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGCCTCAGTTCAGACCTGGGAGGAGTTCAGCCGGGCGGCCGAGAAGCTCTACCTGGCGGACCCCATGA
AGGTACGGGTGGTTCTCAAATACAGGCATGTTGATGGGAATTTGTGTATCAAAGTAACGGATGATCTAGT
TTGTTTGGTGTACAGAACAGACCAAGCGCAAGACGTAAGAAGATTGAGAAATCCACAGTCAGTTAATG
CGACTTATGGTGGCCAAGGAATCCCGCAATGTCACCTATGGAAACAGAA**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_012058
Insert Size: 261 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_012058.3](#), [NP_036188.1](#)

RefSeq Size: 1388 bp

RefSeq ORF: 261 bp

Locus ID: 27058

UniProt ID: [P49962](#)

Cytogenetics: 1 H5

Gene Summary: Signal-recognition-particle assembly has a crucial role in targeting secretory proteins to the rough endoplasmic reticulum membrane. SRP9 together with SRP14 and the Alu portion of the SRP RNA, constitutes the elongation arrest domain of SRP. The complex of SRP9 and SRP14 is required for SRP RNA binding.[UniProtKB/Swiss-Prot Function]