

## Product datasheet for **MC207506**

### Skap2 (NM\_018773) Mouse Untagged Clone

#### Product data:

|                           |   |
|---------------------------|---|
| Product Type:             | Expression Plasmids   |
| Product Name:             | Skap2 (NM_018773) Mouse Untagged Clone  |
| Tag:                      | Tag Free  |
| Symbol:                   | Skap2   |
| Synonyms:                 | 2610021A10Rik; AA960083; BB137539; mSKAP55R; RA70; Saps; Scap2; SKAP-HOM        |
| Mammalian Cell Selection: | Neomycin  |
| Vector:                   | pCMV6-Entry (PS100001)  |
| E. coli Selection:        | Kanamycin (25 ug/mL)  |
| Fully Sequenced ORF:      | >MC207506 representing NM_018773<br>Red=Cloning site Blue=ORF Orange=Stop codon |

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGCCAACCCAGCTGTACCTCTTCTCCCGCCCTCTCCCTGAGGAAATTAGGAACCTGTTGGCAGATG  
TTGAAACATTTGTGCAGACACACTGAAAGGAGAAAATTTATCCAAGAAAGCCAAGGAAAAGAGAGAATC  
TCTCATTAAGAAGATAAAAGATGAAAGTCTGTCTATCTTCAGGAATTTCAAGACAAAGGTGATGCTGAG  
GACGGCGATGAATATGACGATCCCTTTGCTGGCCTGCAGACACGATTTCTTAGCCTCAGAGCGCTACG  
ATAAAGATGATGATGGCCCTCTGATGGAAACCAGTTTCTCCATTGCAGCCCAGGACCTTCCTTTTGT  
CATAAAGGCTGGCTACCTGGAAAAACGCAGAAAAGATCACAGCTTCTGGGGTTTGAATGGCAGAAACGG  
TGGTGTGCGCTCAGCAAAACAGTGTCTATTATTACGGGAGCGATAAAGACAAGCAACAGAAAGGTGAAT  
TTGCAATAGATGGCTATGATGTCAGAATGAACAACACCCTTAGGAAGGATGGAAAGAAAGATTGCTGTTT  
TGAAATCTGTGCTCCCGACAAACGTATCTATCAGTTTACAGCAGCCTCTCCCAAAGATGCTGAGGAATGG  
GTCCAGCAGCTGAAATTTATACTACAAGACCTGGGATCCGACGTTATTCTGAGGATGATGAGGAAAAGAG  
GAGAATTATATGATGACGTTGACCATCTGTGCTGCCGTGACGAGCCCGCAGAGGAGCCAGCCAATCGATGA  
TGAGATTTATGAGGAGCTCCCAGAGGAGGAAGAGGACACTGCTTCAGTGAAGATGGACGACGAAGGGAAG  
GGAAGTCGGACAGTGTGCATCACACCTCAGGAGATAAAAGCACTGATTACGCTAATTTTTACCAGGGTC  
TGTGGGACTGCACTGGAGCTCTGTCTGATGAGTTGTCCTTTAAGCGTGGTATGATTTACATTTTAC  
CAAGGAATACAATAGATATGGCTGGTGGTAGGAGAGATGAAGGGAGCCATTGGCTTGGTGCCTAAAGCC  
TACCTAATGGAGATGATGATATTGA

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

Chromatograms: [https://cdn.origene.com/chromatograms/ja2593\\_e06.zip](https://cdn.origene.com/chromatograms/ja2593_e06.zip)



[View online »](#)

|                               |  |
|-------------------------------|--|
| <b>Restriction Sites:</b>     | Sgfl-Mlul  |
| <b>ACCN:</b>                  | NM_018773  |
| <b>Insert Size:</b>           | 1077 bp  |
| <b>OTI Disclaimer:</b>        | <p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at <a href="mailto:custsupport@origene.com">custsupport@origene.com</a> or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a></p> |
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).   |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>  |
| <b>RefSeq:</b>                | <a href="#">BC003711</a> , <a href="#">AAH03711</a>  |
| <b>RefSeq Size:</b>           | 1656 bp  |
| <b>RefSeq ORF:</b>            | 1077 bp  |
| <b>Locus ID:</b>              | 54353  |
| <b>UniProt ID:</b>            | <a href="#">Q3UNDO</a>   |
| <b>Cytogenetics:</b>          | 6 B3   |
| <b>Gene Summary:</b>          | <p>May be involved in B-cell and macrophage adhesion processes. In B-cells, may act by coupling the B-cell receptor (BCR) to integrin activation. May play a role in src signaling pathway. [UniProtKB/Swiss-Prot Function]</p>  |