

## Product datasheet for **MC207525**

### **Rgs19 (NM\_026446) Mouse Untagged Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Rgs19 (NM\_026446) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Rgs19  
**Synonyms:** 2610042F04Rik; AI324841; AW547781; GAIP  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC207525 representing NM\_026446  
**Red**=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGCCACCCACATGAGGCTGAGAAACAGCACACAGGGCCAGAGGAGGCAGACAGGCCTCCCTCGATG  
CCAGCCATGACGCAGCCCCCTCAGGACCTCCAGTCGCAATCCCTGCTTATGCTGGTGTCTGCTGCTG  
TAGCTGCTCCTGGAACCAAGAACGGCAACGAGCTTGGCAGTTTCTCGGAAAGCAAGCTGCAACCCCTC  
CCCAGCTGTGAAGTCTGCACTCCACCAAGCCCTAAGGAAGTACAGAGCTGGGCACAGTCTTTGACAAGT  
TGATGCATAGCCCCACGGGCCGAGTGTATTCCGGGCATTCTGCGCACAGAATACAGCGAAGAGAACAT  
GCTCTTCTGGCTGGCTTGTGAGGAGTTGAAAGCAGAGGCCAACCAACATGTGGTGGACGAGAAGGCGCGA  
CTTATCTATGAGGACTACGTGTCCATCCTGTCCCCAAGGAGGTAAGCCTGGACTCCCCTGTGCGAGAAG  
GCATCAATAGGAAGATGCAAGAGCCATCGCCACACACATTTGATGATGCACAGCTGCAGATCTACCCCT  
CATGCACCGGGACTCGTACCCTCGATTCCCTTACCTCCCCACCTACCGCTCTCTATTACTCCAGGGGGCC  
CCACAGTCCTCTGAGGCC**TAG**

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_026446  
**Insert Size:** 651 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_026446.4</a></u> , <u><a href="#">NP_080722.1</a></u>
<b>RefSeq Size:</b>	1514 bp
<b>RefSeq ORF:</b>	651 bp
<b>Locus ID:</b>	56470
<b>UniProt ID:</b>	<u><a href="#">Q9CX84</a></u>
<b>Cytogenetics:</b>	2 103.72 cM
<b>Gene Summary:</b>	<p>Inhibits signal transduction by increasing the GTPase activity of G protein alpha subunits thereby driving them into their inactive GDP-bound form. Binds to G-alpha subfamily 1 members, with the order G(i)a3 &gt; G(i)a1 &gt; G(o)a &gt;&gt; G(z)a/G(i)a2. Activity on G(z)-alpha is inhibited by phosphorylation and palmitoylation of the G-protein (By similarity). [UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (2) contains an alternate 5' terminal exon, and it thus differs in the 5' UTR and initiates translation at a downstream in-frame start codon, compared to variant 1. The encoded isoform (b) is shorter at the N-terminus, compared to isoform a. Both variants 2 and 3 encode isoform b.</p>