

## Product datasheet for MC207555

### Sdhc (NM\_025321) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Sdhc (NM\_025321) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Sdhc  
**Synonyms:** 0610010E03Rik; AI316496; AU019277  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC207555 representing NM\_025321  
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGGCTGCGTTCTTGCTGAGACATGTCAGCCGTCCTGCCTCCGAGCCACCTGAATGCTCAGCTTTGTA  
 TCAGAAATGCTGCTCCTTTGGGAACCACAGCTAAGGAGGAGATGGAGCGTTCTGGAAGAAGAACACGAG  
 TTCAAACCGTCTCTGTCTCCCCATTTGACTATCTACAAATGGTCTCTTCTATGGCACTGCTCCGTTTGC  
 CACCGAGGCTCTGGAATAGCCTTGAGTGGAGGGTCTCTCTTTTGGCCTGTCGGCACTGCTGCTTCTG  
 GGAACCTTTGAGTCATATTTGATGTTTGTGAAGTCCCTGTGTTGGGGCCAACACTGATCTACTCGGCTAA  
 GTTTGTGCTTGTCTTCCCGCTCATGTACCACTCACTGAATGGGATCCGACACTTGCTATGGGACCTAGGA  
 AAAGGCCTGGCAATACCCAGGTCTGGCTGTCTGGAGTGGCGGTCTGTTCTTGTCTGTTGTCCTCTG  
 GCGGGCTGGCCGCCCTG**TGA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI

**ACCN:** NM\_025321

**Insert Size:** 510 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_025321.3</a></u> , <u><a href="#">NP_079597.2</a></u>
<b>RefSeq Size:</b>	1179 bp
<b>RefSeq ORF:</b>	510 bp
<b>Locus ID:</b>	66052
<b>UniProt ID:</b>	<u><a href="#">Q9CZB0</a></u>
<b>Cytogenetics:</b>	1 H3
<b>Gene Summary:</b>	Membrane-anchoring subunit of succinate dehydrogenase (SDH) that is involved in complex II of the mitochondrial electron transport chain and is responsible for transferring electrons from succinate to ubiquinone (coenzyme Q).[UniProtKB/Swiss-Prot Function]